USER'S MANUAL

MODEL: LVMCTRL-

LABVIEW MAGNET CONTROL

LabVIEW Driver V2.0

For LabVIEW Ver 8.2 [or later version]

Date Sold:_____

Serial number: _____

PROPRIETARY THIS DOCUMENT CONTAINS CONFIDENTIAL INFORMATION PROPRIETARY TO GMW ASSOCIATES. IT MUST NOT BE REPRODUCED OR DISCLOSED TO OTHERS OR USED IN ANY WAY EXCECPT FOR THE INSTALLATION, OPERATION OR MAINTENANCE OF GMW ASSOCIATES PRODUCTS.

File No: MLVMCTRL GMW Magnet Control Software manual v2.0.doc

Revision Date: Jan 11, 2008

GMW

955 Industrial Road, San Carlos, CA 94070 Tel: (650) 802-8292 Fax: (650) 802-8298 Email: yuqiang@gmw.com Web site: www.gmw.com

Magnet Control S		
Magnet Control S	ystem Basic Requirements	
WARNINGS	[Refer to this section before operation of Electromagnet System]	Section 2
MAGNET CC	NTROL SYSTEM SETUP	Section 3
Hardware	2:	
	Supported field meters	
	Supported power supplies	
	Typical system configurations	
	Digital Teslameter Setup	
	Power Supply Setup Analog Device Interface Setup	
Software:	Analog Device interface Setup	
	Install GPIB-USB driver	
	Install NI-DAQmx software for DAQPad	
	Verify hardware software installation	
	Install magnet control software from GMW	
MAGNET CC	NTROL SYSTEM OPERATION	Section 4
The main	operation screen	
	block diagram	
	nd Stopping the Software	
Operation		
Chart Co		
Log to Fil		
PID Setur		
-	TROL OPERATION	
	troller Overview	
	trol Operation Warning	
	an Open Loop Characterization Plot	
	in Open Loop Characterization Plot	
-	ne PID Controller	
	troller Operation	
	troller Initial Settings Tables	
Open Loo	•	
Setup Sta		
HW Setu		
Wave For		
Table Mo	de.	
MAINTENAN	ICE	Section 5
EXCITATION	I CURVES	Section 6
		2000000
TEST DATA		Section 7
Sorensen	power supply stability test	
Kepco po	wer supply stability test	
APPENDIX		Section 8

 TABLE OF CONTENTS

Section 1

MAGNET CONTROL SYSTEM SPECIFICATIONS

Section 1 MAGNET CONTROL SYSTEM SPECIFICATIONS

Magnet Control System Overview

The Magnet Control System is a special application driver written for using with LabVIEW for Windows version 8.2 or later. This driver provides the user with a GUI [Graphical User Interface] between LabVIEW and the Magnet Control System hardware. It enables operator control of up to two magnet systems comprising Electromagnet, Teslameter and Power Supply in four control modes:

- Open Loop Current Control mode.
- Open Loop Field Control mode.
- Closed Loop Field Control Simple mode.
- Closed Loop Field Control Hybrid mode.

It also supports magnet system without Teslameter in Open Loop Current Control mode or Open Loop Field control mode.

Magnet Control System Basic Requirements

A magnetic field meter from one of the following supported model [Note 1].

- Group3 DTM model: DTM-132/133-DG or DTM-141/151-DG.
 - Senis Magnetic Field Transducer YM/3M [Note 2]
 - MetroLab NMR 2020/2025 with GPIB

A suitable Power Supply which can be either controlled using analog signal [Note 3] or one of the following supported power supplies.

- Danfysik 853/858 power supply with GPIB interface.
- Sorensen SGA with GPIB interface
- Kepco BOP with GPIB interface
- GMW 231HC with NI DAQPad digital control option

A suitable computer system. At the time of printing this manual, suggested computer configuration:

1GB RAM, 1GB available hard drive space, Pentium 4 or equivalent CPU, 2 GHz or above. Operating system Windows XP, SP2 One available USB 2.0 port for GPIB-USB controller if required One available USB 2.0 port for NI DAQPad if required A CD drive for installing software LabVIEW for Windows version 8.2 or later. [Note4]

NI DAQPad 6015/6251.

Required only if using either Senis transducer, GMW 231HC or analog controlled power supply. USB 2.0 interface. Two analog output channels. Output Range: +/-10V. Output Resolution: 16bit. Output impedance: 0.1 Ohm.

Table Mode.

Note:

1. Field (Closed Loop) Control mode cannot be used without a Teslameter installed on the Magnet Control System.

2. Senis Transducer is read through NI DAQ. 3M model is supported but only one field axis is read.

3. Power supply can be bipolar or unipolar. The analog programming range is -10DCV to +10DCV, 16bit resolution.

4. Both source VI and version of stand alone executable program file provided with the CD. The stand alone executable version cannot be altered by the end user. The stand alone executable version does not require LabVIEW.

REFER TO WARNINGS BELOW BEFORE OPERATING ELECTROMAGNET SYSTEM

1 Personnel Safety

If during operation the magnet fringing field is in excess of 0.5mT (5G), it can cause malfunctioning of sensitive electronic and magnetic components. We recommend that warning signs are posted in areas where the fringing field may exceed 0.5mT indicating that a magnetic field may be present.

2 Ferromagnetic Objects

During operation the magnet exerts magnetic attraction towards ferromagnetic objects in the near vicinity of its pole faces. Keep ferromagnetic items clear!

3 Arcing

Magnet stores energy in its field during operation. Do not disconnect any current lead while under load or the magnetic field energy will be discharged across the interruption causing arcing and possible damage to electronic circuits.

4 Coil Hot Resistance

Do not exceed the maximum coil hot resistance for the Electromagnet(s) given in the specifications or coil overheating and possible damage may occur

5 Watches, Credit Cards, and Magnetic Disks

Do not move magnetically sensitive items into the close vicinity of the magnet pole gap. Even some anti-magnetic watches can be damaged when placed in close proximity to the pole gaps during operation. Credit cards, and magnetic disks are affected by magnetic fields as low as 0.5mT (5G). Depending on the previous operating field and the pole gap, the remanent field in the gap can be in excess of 0.5mT (5G) with the magnet power supply off or disconnected.

6 Power off Power Supply

Switch off power supply before performing any hardware setup task around magnet or power supply. The current leads of magnet are isolated from ground. There is an electrical potential difference between the current leads and ground, even when the output current is zero. It could be a safety hazard for personnel or instruments if accidentally contacted to the current leads.

3.1 Hardware.

The Magnet Control System software can be used with different combination of supported field meter and power supply. Some typical combinations of field meter and power supply are listed as below.

-Supported field meter:


Group3 DTM 133/151 with GPIB interface. Senis magnetic field transducer 0.2T/1T/2T/5T with analog output, must be used with NI DAQPad, USB 2.0 interface MetroLab NMR teslameter 2025, GPIB interface.

-Supported power supply:

GMW 231HC, with NI DAQPad, USB 2.0 interface. DANFYSIK 853/858, GPIB/RS232 interface. Kepco BOP, GPIB interface. Sorensen SGA GPIB interface Sorensen SGA analog with NI DAQPad Sorensen SGA analog with GMW Reversal switch, GPIB interface with ADI, USB 2.0 interface with new version of GMW Reversal Switch. Custom power supply with analog programming +/-10V.


Computer interface: USB 2.0 port. GPIB: NI PCI GPIB board NI USB to GPIB controller, using one USB 2.0 port. RS232: USB to RS232 converter.

-Typical system configurations:




System setup 1

Without a field meter, this setup can be used for Open Loop Current and Open Loop Field control mode. Low cost and fast response.




System setup 2

Basic setup used for Closed Loop Control. Low cost, low resolution and can be fast response.




System setup 3

Closed Loop Control setup used for high field stability. Slow response and high resolution.



System setup 4

Two systems control are supported for low speed system setup. Analog control (using NI DAQ) is not supported for two channel control.



System setup 5

3.2 Field meter setup

- 3.2.1 Digital Teslameter [DTM -151/132/133-DG] setup
 - 1. Power down the Digital Teslameter.
 - 2. Remove the bottom cover by undoing the cover securing screw in the center of the cover.
 - 3. The DIP switches located under the bottom cover must be set as detailed in the following table.
 - 4. Replace the bottom cover.
 - 5. Connect the GPIB cable to the GPIB connector located on the rear panel of the Digital Teslameter.
 - 6. Connect the GPIB cable to the GPIB interface board installed in the computer.
 - 7. Connect the Hall Probe.
 - 8. Power up the Digital Teslameter.

Digital Teslameter

(Note 1)

DIP Switch Settings

SW 1-1 OFF	GPIB address
SW 1-2 ON	GPIB address
SW 1-3 OFF	GPIB address
SW 1-4 OFF	GPIB address
	Section 3
SW 1-5 OFF	GPIB address
SW 1-6 OFF	dual primary addressing disabled
SW 1-7 OFF	talker/listener
SW 1-8 OFF	not used
SW 2-1 ON	serviced requests enabled
SW 2-2 ON	EOI asserted
SW 2-3 ON	terminator carriage return
SW 2-4 OFF	double terminator disabled
SW 2-5 OFF	field units tesla
SW 2-6 OFF	no units symbol
SW 2-7 ON	digital filtering ON
SW 2-8 OFF	defaults no action

Note 1: Each DTM requires an unique GPIB address.

3.2.2 MetroLAB NMR 2025 setup

1. Power down the Digital Teslameter.

2. Set Micro-switch on the rear panel.

1 to 5: GPIB address.

6: Set to 0

7: Set to 0. Listener/Talker mode

8: Set to 1. Send <CR><LF> to terminate messages

9: Set to 1. select GPIB (IEE 488) interface.

3. Connect all the cables.

3.3 Power Supply Setup

Below are brief descriptions for setting up powe supply for computer control. For details of how to setup either power supply, please refer to power supply documentation.

3.3.1 GMW 231HC

GMW 231HC has two versions of control. External Analog Control uses a customer provided programming source, such as a Waveform Function Genarator. Internal Analog Control uses a NI DAQ as programming source.

- Ensure the power switch is OFF

- Switch the Control Mode toggle Switch on the back of the GMW 231HC to Internal.

- Connect USB cable from NI DAQPad to computer USB port

- Install all software before power up the DAQPad.

3.3.2 Керсо

- Ensure the power switch is OFF

- Ensure the GPIB Interface board model setup matches power supply model. All new Kepco GPIB boards are set to 200V, 5A model. User needs to change the model on the GPIB board either using a DIP switch (model BIT4882F) or through computer (model BIT4886).

- Set the GPIB address. GMW Control Software uses 1 as default GPIB address. Power supply shipped from GMW will all set GPIB address to 1. If GPIB address needs to be changed, refer to the power supply user manual. Each instrument on the same GPIB bus must have a unique GPIB address.

- Ensure the AC line voltage is correct.

- Connect the currnt cable to magnet. Connect AC cord to power supply.

- Set Kepco power supply front control: Mode to Current, Voltage Control to OFF, Current Control to OFF.

3.3.3 Sorenson

- Ensure the power switch is OFF

- Set the rear panel switch S1-1 to ON to set the unit to remote mode.

- Set the rear panel switch S1-2 to OFF to set the power supply to no SRQ upon power-on.

- Set the rear panel switch S1-3 to ON to set the power supply to be the master device.

- Set the GPIB address. The rear panel switch S1-4 to S1-8 is GPIB address. Power supply shipped from GMW will all set GPIB address to 1. If GPIB address needs to be changed, refer to the power supply user manual. Each

instrument on the same GPIB bus must have a unique GPIB address.

- Ensure the AC line voltage is correct.

- Connect the currnt cable to magnet. Connect AC cord to power supply.

3.3.4 Danfysik

- Ensure the power switch is OFF
- Set the GPIB address. If using NI GPIB-232CV-A, the following needs to be check/set.
- 1. Danfysik power supply main control board:
 - Remote Line needs to be set to RS232: Jumper on ST9. ST10 and ST11 must be left open without jumper.
 Remote Line communication parameters are as Danfysik power supply defult: 9600 BAUD, No Parity, 8 Data bit, 1 Stop bit.
- 2. Cable between main control board and GPIB-232CV-A must be null modem cable, DB25 Male/DB9 Female.
- 3. NI GPIB-232CV-A.
 - Mode: D
 - Config:SRQ Disabled, small serial buffer
 - GPIB address: 1
 - Hand Shake:XON/OFF disabled
 - Data formate: 8 Data bit/No Parity/1 Stop bit
 - Baud: 9600

Power supply shipped from GMW will all set GPIB address to 1. If GPIB address needs to be changed, refer to the power supply user manual. Each instrument on the same GPIB bus must have a unique GPIB address.

- Ensure the AC line voltage is correct.
- Connect the currnt cable to magnet. Connect AC cord to power supply.
- Power on the Danfysik power supply. From front panel, change the Local/Remote control to REMote.

3.4 Analog Device Interface Setup for NI DAQPad [6015]

- 1. Ensure the power switch is OFF
- 2. Connect NI DAQPad analog output (BNC, male) to power supply analog control input.
- 3. Connect USB cable from NI DAQPad to computer USB port.
- 4. Install all software before power up the DAQPad.

3.5 Software Installation 3.5.1 Install GPIB-USB driver

- 1. Insert NI-488.2 CD from National Instruments into CD-ROM drive.
- Select "Install Software" from pop-up window.
 ³⁷ NI-488.2 for Windows



3. Click on "Next".

NI-488.2 for Windows		
Welcome to NI-488.2 for Windo		
It is strongly recommended that you exit all programs b	before running setup.	
This program is subject to accompanying National Ins License Agreement.	This program is subject to accompanying National Instruments Software License Agreement	
For patents covering National Instruments products, refer to the appropriate location: Helps >Patents in your software, the patents.txt file on your CD, or ni.com/patents.		
NI-488.2 [™]		
Version 2.4		
ni.com/gpib	©2005 National Instruments. All rights reserved.	
	<< Back Next >> Cancel	

4. Select destination directory. Use default setting recommend. Click on "Next".

3.2 for Windows		
Destination Directory Select the primary installation	directory.	
Instruments\NI-488.2 folder. To	are will be installed by default to th rstall into a different folder, click th oses individual component location next dialog.	he Browse button,
Instruments\NI-488.2 folder. To and select another folder. To ch	nstall into a different folder, click th pose individual component location	he Browse button,

5. Select Installation option. <u>Select "Custom". Click on "Next".</u>

ANI-488.2 for Windows	
Select Installation Option Select one of the following options.	
Typical Use this option to install the most common application features Use this option to install the most common application features Use this option to select the application features you want to advanced users. NI-NISA Software is required if you want to run instrument or choose this option in order to install NI-VISA.	to install. Recommended for driver software. You must
<< <u>B</u> ack	<u>N</u> ext >> <u>C</u> ancel

6. Select Features to install. Click on "NI-VISA 3.2". Select "Install this feature to a local drive" from the drop-down menu. Click on "Next".

🗏 NI-488.2 for Windows	
Features Select the features to install.	
NI-468.2.2.4 Driver and Applications GPIB Analyzer Language Support Microsoft C/C++ Borland C/C++ Visual Basic LabV/EW LabV/EW LabV/EW LabV/Indows/CVI Other NIT Framework 1.1 Langua	National Instruments VISA driver version 3.2. VISA provides a API for controlling V/3, GPIB, Serial, P/3 and other types of instruments.
Tinstall this feature to a local drive for 3.1.	This feature and its selected subcomponents will require 7.58 MB of disk space.
Restore Defaults	Cancel

7. Accept the License Agreements. Click on "Next".


NI-488.2 for Windows	
License Agreement You must accept the license(s) displayed below to proceed.	
NATIONAL INSTRUMENTS SOFTWARE LICENSE AGREE	MENT
INSTALLATION NOTICE: THIS IS A CONTRACT. BEFORE YOU DOWNL SOFTWARE AND/OR COMPLETE THE INSTALLATION PROCESS, CAI THIS AGREEMENT. BY DOWNLOADING THE SOFTWARE AND/OR CL APPLICABLE BUTTON TO COMPLETE THE INSTALLATION PROCESS TO THE TERMS OF THIS AGREEMENT AND YOU AGREE TO BE BOU AGREEMENT. IF YOU DO NOT WISH TO BECOME A PARTY TO THIS AND BE BOUND BY ALL OF ITS TERMS AND CONDITIONS, CLICK TH BUTTON TO CANCEL THE INSTALLATION PROCESS, DO NOT INSTAL SOFTWARE, AND RETURN THE SOFTWARE WITHIN THIRTY (30) DAN OF THE SOFTWARE (INCLUDING ALL ACCOMPANYING WRITTEN MA WITH THEIR CONTAINERS) TO THE PLACE YOU OBTAINED THEM. AI SHALL BE SUBJECT TO NI'S THEN CURRENT RETURN POLICY.	REFULLY READ JCKING THE , YOU CONSENT IND BY THIS AGREEMENT E APPROPRIATE L OR USE THE /S OF RECEIPT .TERIALS, ALONG
 I do not accept the Licen 	se Agreement(s). e License Agreement(s).
(<u>Rack</u>)	>> <u>C</u> ancel

8. Start Installation. Click on "Next".


NI-488.2 for Windows	
Start Installation Review the following summary before continuing.	
Features to be Added or Changed •NI-488.2.2.4 Driver and Applications GPIB Analyzer Language Support Microsoft C/C++ Borland C/C++ Visual Basic LabVIEW LabWindows/CVI Other • NI Measurement & Automation Explorer 3.1.1	
Click the Next button to begin installation. Click the Back button to change the insta	allation settings.
(< <u>B</u> ack)	> <u>C</u> ancel

9. Wait until the installation complete. Click on "Next".

邊 NI-488.2 for Windows	
Overall Progress	
Currently installing NI-488.2 2.40. Part 2 of 15.	
	<< Back Next >> Cancel



10. Select GPIB device: GPIB-USB-HS. Click on "Next"




11. Select "Restart. Click on "Finish".

🕅 Add GPIB Hardward	e Wizard 📃 🗖 🔀
NI-488.2 [™] High Performance GPIB Software	The GPIB interface has been added to your system. For additional help, refer to the GPIB Hardware Guide. <u>View GPIB Hardware Guide</u> To add additional GPIB interfaces, click Add Another Interface. To install a plug-in board (e.g. PCI-GPIB), select Shutdown and click Finish. To install an external interface (e.g. GPIB-ENET/100) or if the plug-in board is already in the computer, select Restart and click Finish. • Shutdown • Estart
	Add Another Interface Finish

12. After restart computer, connect GPIB-USB controller to computer USB port. Found New Hardware Wizard window appears. Select "No, not this time." Click on "Next"



13. Select "Install software automatically". Click on "Next"



14. Select "Search for the best driver in these locations". Click on "Next".

Found New Hardware Wizard		
Please choose your search and installation options.		
Search for the best driver in these locations.		
Use the check boxes below to limit or expand the default search, which includes local paths and removable media. The best driver found will be installed.		
Search removable media (floppy, CD-ROM)		
Include this location in the search:		
D:\ Browse		
O Don't search. I will choose the driver to install.		
Choose this option to select the device driver from a list. Windows does not guarantee that the driver you choose will be the best match for your hardware.		
< <u>B</u> ack <u>N</u> ext > Cancel		

15. Wait until finished. Click on "Finish".

Found New Hardware Wizard		
	Completing the Found New Hardware Wizard	
	The wizard has finished installing the software for:	
	The wizard has finished installing the software for:	
	Click Finish to close the wizard.	
	< Back Finish Cancel	

16. Select the option for registration hardware. Click on "Finish".


National Instrument	ts Registration Wizard			
	The National Instruments Registration Wizard helps you discover and register your National Intruments products to facilitate technical support and ensure you receive important updates.			
	red the following National Instr discovered, you can add them		cts in your system. If you have	
	Product Name	Serial Number	Refresh	
	GPIB-USB-HS	0120f558	Add More Remove	
Would you like to regi	ster these products?			
	🔘 Register Now			
	🔵 Remind me again in 14 d	lays		
	💿 Do not register any of m	ny products		
How does National In	struments use registration info	ormation?		
		<< Back	Finish Cancel	

17. "NI 488.2 Getting Started Wizard" window appears. After verifying hardware and software installation, click on "Exit".

🕅 NI-488.2 Getting Started Wizard			
Click on the following steps to use the NI-488.2 Software for Windows with your GPIB instruments:			
Verify your hardware and software installation			
Communicate with your instrument			
✓ Do not show at Windows startup	Exit		

3.5.2 Install NI-DAQmx software for DAQPad.

Following steps are for NI-DAQmx 8.0.1. Actual step and screen may vary depend on version of DAQmx 1. Insert DAQmx CD from National Instruments. From the pop-up window, select "Install NI-DAQmx/VI Logger".



2. Product Information window appears. Click on "Next".5

鴉	NI-DAQ 8.0.1
	Product Information Please read the following information about the products to be installed.
	 NI-DAQ 8.0 Device Support NI-DAQ 8.0 supports the latest DAQ devices, including M Series, E Series, S Series, B Series, USB Series, the NI PCI-7041/80400E for LabVIEW RT, Analog Output/Waveform, Counter/Timer, Digital 1/0, DSA, SCAL, and SCC devices. Refer to the NI-DAQ Readme for supported device models. Traditional NI-DAQ (Legacy) is an older driver with an outdated API that supports some additional older DAQ devices. ApI to rot all of the latest DAQ devices. ApI closed with earlier versions of NI-DAQ will run under Traditional NI-DAQ (Legacy). You can use both NI-DAQmx and Traditional NI-DAQ (Legacy) in the same computer, but must reset devices before changing from Traditional NI-DAQ (Legacy) to NI-DAQ will run NI-DAQ Readme for guidelines. Support Removed from NI-DAQ 7.0 and later NI-DAQ 7.0 and later supports different devices, operating systems, application software, and programming language versions than previous releases of NI-DAQ. To Unity DAQ, the DAQC and 500/700 family, parallel port DAQPads, and several other older devices. Do not install NI-DAQ 7.0 or later if your existing applications include unsupported components. Refer to the NI-DAQ Readme file on the distribution.
	Save File << Back Next >> Cancel

3. Select destination directory. Recommend using default directory setting. Click on "Next"

掲 NI-DAQ 8.0.1	
Destination Directory Select the primary installation directory.	NATIONAL INSTRUMENTS
Destination Directory C:\Program Files\National Instruments\	Browse
	<< Back Next >> Cancel

4. Select features to install. Do not change default settings. Click on "Next".

NI-DAQ 8.0.1	
Features Select the features to install.	
NI-DAQmx 8.0 LabVIEW 8.0 Support LabVIEW 7.1 Support LabVIEW 7.1 Support LabVIEW 7.0 Support LabVIEW 7.0 Support LabVIEW 7.0 Support LabView 7.0 Support NET Framework 1.1 Languages Su Measurement Studio for Visual C++ 2 Microsoft Visual Basic 6.0 Support VI Logger VI Logger NI Measurement & Automation Explorer 4.0	NI-DAQ mx, the next-generation NI-DAQ driver with new VIs/functions and development tools for programming your measurement devices, has new features and advantages over the previous version. $\overbrace{VVV} \bigoplus \bigoplus \bigoplus \bigoplus$ This feature will be installed on the local hard drive. This feature and its selected subcomponents will require 137 MB of disk space.
< >>	
Directory for NI-DAQmx 8.0	
C:\Program Files\National Instruments\NI-DAQ\	Browse
Restore Defaults Disk Cost	<pre></pre>


5. Accept the License Agreements. Click on "Next".

₩ NI-DAQ 8.0.1	
License Agreement You must accept the license(s) displayed below to proceed.	
NATIONAL INSTRUMENTS SOFTWARE LICENSE AGREEM	ENT 🔼
INSTALLATION NOTICE: THIS IS A CONTRACT. BEFORE YOU DOWNLO SOFTWARE AND/OR COMPLETE THE INSTALLATION PROCESS, CAR THIS AGREEMENT. BY DOWNLOADING THE SOFTWARE AND/OR CLU APPLICABLE BUTTON TO COMPLETE THE INSTALLATION PROCESS, TO THE TERMS OF THIS AGREEMENT AND YOU AGREE TO BE BOUN AGREEMENT. IF YOU DO NOT WISH TO BECOME A PARTY TO THIS A AND BE BOUND BY ALL OF ITS TERMS AND CONDITIONS, CLICK THE BUTTON TO CANCEL THE INSTALLATION PROCESS, DO NOT INSTALL SOFTWARE, AND RETURN THE SOFTWARE WITHIN THIRTY (30) DAYS OF THE SOFTWARE (INCLUDING ALL ACCOMPANYING WRITTEN MAT WITH THEIR CONTAINERS) TO THE PLACE YOU OBTAINED THEM. ALL SHALL BE SUBJECT TO NI'S THEN CURRENT RETURN POLICY.	EFULLY READ KING THE YOU CONSENT D BY THIS GREEMENT APPROPRIATE OR USE THE OF RECEIPT ERIALS, ALONG
1. <u>Definitions.</u> As used in this Agreement, the following terms have the fol	
A "You " Means you the individual using the SOFTWARE as well as accept the License	
OI do not accept the L	
<< Back Next >>	

6. Start Installation.



7. Wait until finished. Click on Next". Click on "Restart" when prompt.



- 8. After restart computer, connect DAQPad to computer USB port.
- 9. "Found New Hardware Wizard" window appears. Select "No, not this time." Click on "Next".

Found New Hardware Wizard	
Hard Viologing the Win	come to the Found New ware Wizard will search for current and updated software by on your computer, on the hardware installation CD, or on dows Update Web site (with your permission). <u>ur privacy policy</u>
Softward	ndows connect to Windows Update to search for ? es, this time only ies, now and gvery time I connect a device to, not this <u>time</u> ext to continue.
	< <u>Back</u> Next > Cancel

10. Select "Install the software automatically". Click on "Next".

Found New Hardware Wiz	ard
	This wizard helps you install software for: DAQPad-6015 If your hardware came with an installation CD or floppy disk, insert it now. What do you want the wizard to do? Install the software automatically (Recommended) Install from a list or specific location (Advanced)
	Click Next to continue.
	< <u>B</u> ack <u>N</u> ext > Cancel

- 11. Wait until finished. Click on "Finish".
- 12. NI-Max will detect new data acquisition device. Select "Take No Action". Click on "OK".

🕅 New Data Acquisition Device	×
NI-DAQmx has detected a new data acquisition device:	
🏋 DAQPad-6015	
What would you like to do?	
Run Test Panels	
Configure This Device Using MAX	
Run VI Logger	
Take No Action	
🔲 Always Take This Action	
OK Cance	

3.5.3 Verify hardware software installation

- 1. Open "Measurement and Automation Explorer" from National Instruments.
- 2. Go to "My System/Device and Interfaces". There should be two entries, one as "GPIB0(GPIB-USB-HS)", one as "NI-DAQmx Devices"



3. Right click on "GPIB0 (GPIB-USB-HS)". Select "Scan for Instruments" from the drop down menu. All GPIB devices on the GPIB bus should appear in the right side window. The GPIB addresses will be needed later for software setup.

e Edit. Yiew Tools Help			
Configuration 🔗	Properties Scan Fo	or Instruments 🛛 😭 NI Spy 🧉	& GPIB Analyzer
My System	Name	Value	Description
	FGPIB Interface Number	1 Yabe 0 Primary Address = 1 Primary Address = 2	Description

4. Click on "NI-DAQmx Devices". The DAQPad should show on the right side window, without red cross. Note the device reference is "Dev1". If more than one DAQPad have been connected to this computer, there will be more devices showing. This reference for device will be needed later for software setup.

💐 NI-DAQmx Devices - Measurement	A Automation Explorer		
	a Automation Explorer		
<u>File Edit View Tools H</u> elp			
Configuration			
Signed System Image: System Image: System Image: System	NI-DAQmx Devices ₩ NI DAQPad-6015: "Dev1"	Device Type	

3.5.4 Install magnet control software from GMW.

- 1. Insert software CD from GMW. Select D:\Installer\Setup.exe to start the installation.
- 2. Select destination directory. Make changes if wanted. It is recommended leave the "Directory for National Instruments products" unchanged. Click on "Next"

😼 Savannah Copley Power Supply Control ver 1.0
Destination Directory Select the primary installation directory.
All software will be installed in the following location(s). To install software into a different location(s), click the Browse button and select another directory.
Directory for Savannah Copley Power Supply Control ver 1.0 C:\Program Files\Savannah Copley Power Supply Control ver 1.0\ Browse
Directory for National Instruments products C:\Program Files\National Instruments\ Browse
<< <u>B</u> ack Next >> Cancel

3. Accept the License Agreements. Click on "Next"

📲 Savannah Copley Power Supply Control ver 1.0
License Agreement You must accept the license(s) displayed below to proceed.
NATIONAL INSTRUMENTS SOFTWARE LICENSE AGREEMENT
INSTALLATION NOTICE: THIS IS A CONTRACT. BEFORE YOU DOWNLOAD THE SOFTWARE AND/OR COMPLETE THE INSTALLATION PROCESS, CAREFULLY READ THIS AGREEMENT. BY DOWNLOADING THE SOFTWARE AND/OR CLICKING THE APPLICABLE BUTTON TO COMPLETE THE INSTALLATION PROCESS, YOU CONSENT TO THE TERMS OF THIS AGREEMENT AND YOU AGREE TO BE BOUND BY THIS AGREEMENT. IF YOU DO NOT WISH TO BECOME A PARTY TO THIS AGREEMENT AND BE BOUND BY ALL OF ITS TERMS AND CONDITIONS, CLICK THE APPROPRIATE BUTTON TO CANCEL THE INSTALLATION PROCESS, DO NOT INSTALL OR USE THE SOFTWARE, AND RETURN THE SOFTWARE WITHIN THIRTY (30) DAYS OF RECEIPT OF THE SOFTWARE (INCLUDING ALL ACCOMPANYING WRITTEN MATERIALS, ALONG WITH THEIR CONTAINERS) TO THE PLACE YOU OBTAINED THEM. ALL RETURNS SHALL BE SUBJECT TO NI'S THEN CURRENT RETURN POLICY. 1. Definitions. As used in this Agreement, the following terms have the following meanings:
 I accept the License Agreement(s). I do not accept the License Agreement(s).
<< Back Next >> Cancel

4. Start Installation. Click on "Next".

🖁 Savannah Copley Power Supply Control ver 1.0	
Start Installation Review the following summary before continuing.	
 Upgrading National Instruments system components Adding or Changing Savannah Copley Power Supply Control ver 1.0 Files 	
, Click the Next button to begin installation. Click the Back button to change the installation settings.	
Save File) << Back	ncel

5. Wait until finished. Click on "Finish".

😴 Savannah Copley Power Supply Control ver 1.0	
Installation Complete	
The installer has finished updating your system.	
	Finish

6. Restart computer when prompted.

4.1 The main operation screen

- Shows Field readings from the field meter if connected and configured correctly.

- Window control for selecting Window Sub Panels, Operation, PID, Chart, Wave Form, Table Mode and Ramp Mode

- Shows the system status.

- Access to Hardware setup menu. (HW Setup button)


- Controller Mode selects from one out of four control modes: Open Loop Current, Open Loop Field, Closed Loop Field Simple or Closed Loop Field Hybrid.

- STOP button programs the power supply to zero current and stops the program exection.

- System Parameters Strip chart shows system parameters on a strip chart under Chart tab. More details later.

GMWAssociates		PID	Chart	Wave Form	Table M	ode Ramp	Mode
Software Version: 2.02 Nov 14, 2007	Operation Mode Open Loop Current Control Open Loop Field Control	bl					
Esc) STOP	Close Loop Field Control S						-
HW Setup	Field Probe Location						
Log to File Remote Mode	OProbe at Secondary po						
Log to File:	Field Meter Range			Set Field	÷) C	.00000) (T)
Remark: Your remark here.	0.01T	г		Set Curre	ent ⁽)	0.0000) (A)
Start Log Manual Log on Manual Log			Pro	obe Field	, Вр 🚺	.00000	(T)
Sample Interval	Custom		Sam	nple Field	, Bs 【	.00000	(T)
System Status			Οι	itput Cur	rent	0.00	(A)
	To Be Add		-100	-80 -60 -40	-20 0 2	0 40 60	80 100
				tput Volt		0.00	(V)
Loops Read loop 0			20 17.				

4.2. Software block diagram.



Software block diagram

4.3. Starting and Stopping the Software

To open the program, go from Start>>Program>GMW Magnet System Closed Loop Control>GMW Magnet System Closed Loop Control V2.

To start the software, click on the arrow button on the top of the screen. To stop the software, click on Stop button or use <Esc> key.

The software will load default settings from a configuration file (Config.txt) every time the software starts. To change the default settings click on HW Setup.

For first time setup, click HW Setup before starting the software. It will open the hardware setup window.

If more than one power supply/field meter are enabled in hardware setup, another tab will appear to show the second system. The second system display and control is identical with the first system. The description below is for the first system unless specified.

4.4. Operation Menu.

4.4.1 If controlling more than one power supply, power supply under control is selected by click on the proper tab. If no second magnet tab shown, the software is controlling one power supply.

4.4.2 Operation Mode:

Open Loop Current Control. Controls power supply current output by setting the Set Current to desired current value. User can not change Set Field under Current mode.

To change the current output of power supply under current mode, select Open Loop Current Mode, input desired current at Set Current. Output current will change after hit <Enter> or click anywhere outside the Set Current box.

Open Loop Field Control. Set power supply current output to reach a given field. Software will determine the output current from the given field through a predefined Field vs. Current excitation curve, (OLC, Open Loop Characterization file). The OLC file can be made either using a field meter support or a simulation file. OLC file in use is displayed under PID tab. The accuracy of Open Loop Field Control depends on the accuracy of the OLC file.

To change the current output of power supply under open loop field mode, select Open Loop Field Control, input desired field at Set Field. Output current will change after hit <Enter> or click on any where outside the Set Field box.

Close Loop Field Control Simple. Adjust power supply current output to maintain a given field by setting the Set Field to desired field value. Software will compare the desire field and the reading of field meter, adjust the output current accordingly. User can not change Set Current under Field mode.

To change the current output of power supply under close loop control field simple mode, select Close Loop Field Control Simple, input desired field at Set Field. Output current will change after hit <Enter> or click any where outside the Set Field box.

Close Loop Field Control Hybrid. In this mode the controller uses data from the Open Loop Characterization File and operates in PID mode only when the controller is within a 5% window of the desired Field Setpoint. The Open Loop Characterization File provides the controller with the field versus current excitation curve. The controller will not operate correctly if an Open Loop Characterization file is not loaded or the file has incorrect data for the magnet. When the user gives the system a new Field Setpoint the controller calculates an initial current setting for the desired Field Setpoint and transmits it to the Power Supply. During the ramp up or down to the desired Field Setpoint the PID Controller is disabled. When the actual field reaches a value within the Hybrid control window (defined in the HW setup and loaded when software started) the PID controller is activated. During the change over to PID control the Derivative control is disabled to allow the controller to settle before taking control of the system.

Generally the Hybrid PID mode is preferred since it provides the fastest transition time to a new Field Setpoint. User can not change Set Current under Field mode.

To change the current output of power supply under close loop field control mode, select Close Loop Field Control Hybrid, input desired field at Set Field. Output current will change after hit <Enter> or click on any where outside the Set Field box.

4.4.3 Field Probe Location:

Probe at Sample Position. Field reading from probe is used directly to control the field without change. The software is controlling the field at the point the probe is located.

Probe at Secondary Position. Field reading from probe is modified before used for controlling the field. The modification is made either by a constant multiplier or looking into two OLC files at Sample position and Secondary position. This allows controlling field indirectly if the probe cannot be placed at the intend controlling location. The constant multiplier and the two OLC files are defined in HW Setup window.

4.4.4 Closed Loop Field Control Operation Warning

Before attempted to run the Magnet Control System in the Closed Loop Field control mode check the items listed below.

- 1 The Teslameter is powered up.
- 2 The Teslameter probe is inserted in the magnet gap securely.
- 3 The Teslameter is set to the correct field range.
- 4 The Teslameter is connected to the GPIB bus and is returning correct field readings.
- 5. The Teslameter is reading the field correctly.

Running close loop field control without the field meter properly installed will cause the power supply to ramp up to max current. It is dangerous for both the operator and equipment setup.

4.4.5. Field Meter Range. Select the field meter full-scale range to be higher than the maximum possible field.

4.4.6. Field Reading.

- Probe Field. Shows the field meter reading.

- Sample Field. Shows the field at sample position. If the probe is located at sample position, the Probe and Field and Sample Field are the same. If the probe is located at secondary location, the Sample Field is calculated field.

4.5. Chart Control.

- Auto Scale ON/OFF: Turn on/off auto scaling for Y-axis.
- Time Axis: Set the time span for X-axis.

- Channel Select: Select/deselect which trace to be shown on the chart. This setting also determines if the trace will be saved to a file.

- **4.6.** Log to File. Save data to a text file. The traces to be saved are selected from the Chart display, Channel Select. If a trace is unselected, that trace will not be saved to file. After Log to File is started, Select/unselect trace will cause software misplace data. User should avoid changing Change Select during logging.
 - Log To File: Shows the data file name/path.
 - Start Log: Click on Start Log. A window will prompt user for file name/directory to be saved.
 - Auto/Manual: If Auto is selected, data will be saved every time interval defined by Sample Interval. If Manual is selected, data will be saved every time the user clicks on Manual Log.

4.7. PID Setup.

This menu is used for tuning the PID Controller

Before using the PID controller for the first time the controller will need to be tuned to obtain optimum performance. Follow the instructions listed below.

- 1. Select the PID Setup tab. The PID control panel will appear.
- 2. Set the Controller Mode required, either Hybrid PID or Simple PID control.
- 3. Enable the display of the field on the stripchart by selecting it on the Channel Select on the side of the stripchart.
- 4. Set the PID and Lead Lag controls to the initial values. Set the Lead Lag gain to 1, lag time to 0.2, lead time to 0.
- 5. Check that the correct Open Loop Characterization file is loaded. If no Open Loop Characterization file is found, a Default.olc will be loaded. [See Section 4.10 for Loading an Open Loop Characterization Plot.]
- 6. Set the Integral and Derivative controls to 0.00.

7. Enter a Field Setpoint value of 0.100000 Tesla. Retry with a Field Setpoint of 0.200000Tesla.

8. Increase lag time until the field can reach the target field quickly and does not have overshot.

- 9. Observe the trace of the field for the controller response on the System Parameter chart. It should not oscillate.
- 10. Increase the PID Gain a small increment.
- 11. Repeat steps 9 and step 10.
- 12. Continue repeating steps 10 and 11 until the controller oscillates in a steady state. Steady state oscillations mean the oscillations neither decrease or increase over a period of time.
- 13. Reduce the PID Gain a small increment, repeat steps 8 and step 9. If the controller does not oscillate then the system is set at maximum gain.
- 14. Reset the PID Integral and Derivative controls to the initial values given in Section 9 of this manual.
- 15. Repeat step 8 and step 9. If the system oscillates then decrease both the Integral and Derivative controls. If the systems overshoot the Field Setpoint then increase both the Integral and Derivative controls.

Note: 1. Integral and Derivative controls should have an approximate ratio of 4 to 1 between their settings.

[e.g. If the Integral control is set to 6.00 then the Derivative should be set to 1.50].

2. All the PID tuning control settings are saved only at HW Setup.

4.8 Open Loop File.

4.8.1 Load Open Loop File

Before the PID controller is run, a B v I Open Loop Characterization File must first be loaded that matches the magnet and gap being used. This B v I Open Loop Characterization file (OLC file) can be generated either using a field meter or using a simulated field vs. current curves. An example OLC file below shows the file format. Current and Field columns is TAB delimited.


Current (Amps)	Field (Tesla)
0.0000E+0	-3.0350E-4
8.3333E-1	6.4058E-2
1.6667E+0	1.2882E-1
2.5000E+0	1.9346E-1
3.3333E+0	2.5757E-1
4.1667E+0	3.2085E-1
5.0000E+0	3.8263E-1


If the appropriate OLC file exists in the system, the OLC file can be loaded as detailed below. If the appropriate file does not exist the user will need to make the Open Loop Characterization File, [see Section 4.8.2 "Making an Open Loop Characterization Plot"]. To make the OLC to be loaded every time the software starts, making the OLC file in use as default file in the HW Setup window. [see Section 4.10 "HW Setup"]

1. Click on the LabVIEW Run arrow to start the Magnet Control software [if not already running].

2. Click on the PID tab and then click on Open Loop Characterization.

CMIU.		
GMWAssociates Magnet System Closed Loop Control Software Version: 2.0a Oct 27, 2007	Operation Menu PIE	Chart Ramp Mode Table Mode Setup Status
(Esc) STOP HW Setup	lag time (mm) (0.005 lead time (mm) (0.001	entegral time (Tri, mor) § 0.003 derivative time (Tri, mor) § 0.003
Log to File Remote Mode	-Open Loop File in Use: Secondar	at Centrel 1/2.a/1470-40-20_1-1-2007.de
StartLog J Manual Log on J Manual Log Sample Interval 1 (sec)	Kepco	Probe Field 0.00000 (T) Sample Field 0.00000 (T)
System Status		Output Current 0.00 (A)
		Output Voltage 0.00 (V)





- 4. Select an Open Loop Characterization file by clicking on the desired file name. [A valid Open Loop Characterization file will have .olc extension on the end of the filename].
- 5. Click on OK or press the Enter key. The selected file will be loaded.
- 6. Check the selected file name is shown in the Loaded Open Loop Characterization File indicator. If correct then click on the Done button. The Open Loop Characterization Plot panel will close.
- 7. To use the selected OLC file each time the software starts, input the OLC file name and path in HW Setup>Power Supply>Open Loop File Path. Save the configuration.


4.8.2 Making an Open Loop Characterization Plot.

Follow the previous steps to open the Open Loop Characterization screen.

- 1. Set the Max and Min current for the intend power supply operating range. It is not necessary the Max and Min current match the output range of the power supply.
- 2. Set the Steps. More steps will generate more accurate Open Loop Characterization Plot.
- 3. Click on Begin button.
- 4. The system will start taking data and report status in the status indicator. When the system reports it is Idle then click on the Save Button. Name the plot file using the correct format in the Filename box and click on OK or press the Enter key. [See next Section 4.8.3 for description of Open Loop Characterization filename format].

4.8.3 Open Loop Characterization filename format

It is recommended the Open Loop Characterization File name has the following format.



4.9. Setup Status.


4.10. HW Setup.

HW Setup is used for configurating system hardware. All setup parameters are saved to a configuration file [Config.txt]. Software will load the settings from Config.txt every time the software is started. The user must stop and restart the software after the change of settings to enable the new settings to be applied.

4.10.1. Click on the HW Setup button. The HW Setup menu will appear.

Save Configuration	GMW Magnet Con Date Time	Software Title	Software Version 2.0a
Load Configuration Exit Setup	10/29/2007 2:17:53 PM	GMW Magnet Control Software Config file path H://uqiang's Data/Lab/IEW code/GMW Magnet Control V2.a/config.txt	Software Release Date 10/14/2007 Configuration File Date/Time
tup Yower Supply 1 Power Supply 1 Power Supply 1 Setup	2 DAQPad Setup System Se	etup PID Control Setup 1	
YES Power Supply Tockpl YES Power Supply Type 2231 HC Power Supply Voltage Range (V) \$150 Min Current Range (A) \$20 Max Current Range (A) \$20 Power Supply Interface \$UISB DAQ DAQmx Output Channel \$20 Yes	No Field Meter Enabled? Select Field Meter Brand: Senis Probe Max Range I T Field Meter inteface: DAQmx Input Channel DAQmx Input Channel	PID pains proportional gain (Kc) \$ 5.000 integral time (Tr, min) \$ 0.011 derivative time (Td, min) \$ 0.001 output range \$ 20.000 output tow \$ -20.00 Lead-Lag tuning parameters gain \$ 5.000 lag time (min) \$ 0.001 lead time (min) \$ 0.001	Primary Open Loop Fie path 1 C: Program Files/GMW Magnet Control V2.4(Test OLC 5201 Kepco Sensi J002207.dc Use single Inear factor (F-Bsample/Bprobe) Use Secondary Open Loop fie Probe Scale Factor 1.23

- 4.10.2. Power supply 1 and power supply 2 have identical setup panels.
- 4.10.3. Select Power supply Type.
- 4.10.4. Input the minimum and maximum current range according to the power supply.
- 4.10.5. Select Power Supply Interface.
- 4.10.6. Select DAQmx Output Channel from the pull-down menu if using 231 HC or customer power supply with analog control. This is the analog output channel from National Instruments DAQPad. The Dev is the device address automatically assigned by software. DAQPad has two analog output channels, AO 0 and AO 1. Set the channel according to the physical connection between DAQPad and power supply. If only one power supply is under control, the default output channel is AO 0. If no pull-down menu is available, check the NI software installation and verify the hardware according to the Software/Hardware Installation section.
- 4.10.7. Select Field Meter Brand.
- 4.10.8. Select Field Meter Model or Range.
- 4.10.9. Select Field Meter interface.
- 4.10.10. Select the Field Meter GPIB addresses or DAQPad Channel from the pull-down menu. All devices on the same GPIB bus must have an unique GPIB address. The pull-down menu will show all the devices on the same GPIB interface. If no device shows on the pull-down menu, check the GPIB cable connection and power connection to the devices.
- 4.10.11. Input the PID settings. If first time setup, leave the PID settings unchanged.
- 4.10.12. Click on the button beside the Open Loop File path. Select the correct Open Loop File from the pop-up window.
- 4.10.13. On DAQPad setup, ensure the output sample rate is below 30.
- 4.10.14. Click on the System Setup tab.



- Read Field average number: Samples used for calculating moving average of the field meter reading.
- Max Output Slew Rate: Sets the limit for maximum current ramp rate.
- Hybrid Window: Window used for Hybrid PID control.
- 4.10.15. Save Configuration. Click on Save Configuration button. All the settings will be saved to file Config.txt. Each time Save Configuration is clicked, old settings are overwritten by the new settings. To preserve the old settings for latter use, copy the Config.txt file and rename the file using Window Explorer. The software has to be restarted to apply the new settings.

4.10.16. Load Configuration. Click on Load Configuration will let user load from a previous saved configuration file other than default Config.txt.

4.11. Wave Form Mode.

Wave Form mode will generate a wave form according to user settings and use the wave form to set the output current. Operate on current mode only.

4.11.1. Select Ramp Mode tab.

GMWAssociates Magnet System Closed Loop Control	Operation Menu P	ID Chart Wave Form Table Mode Ramp Mode
Software Version: 2.02 Nov 14, 2007 (Esc) STOP		Ramp Waveform ON \$Guare Wave OFF Offset(A) \$0
Log to File Remote Mode		Frequency(Hz) \$0.5 Amplitude(A) \$0 Single Ramp
Remark: Your remark here.		Probe Field, Bp 0.00000 (T)
Log on Sample Interval	Custom	Sample Field, Bs 0.00000 (T)
System Status No error	To Be Add	Output Current 0.00 (A)
Loops		Output Voltage 0.00 (V)

4.11.2. Select desired waveform.

hart	Ramp Mode	Table Mode	Set
	Ramp ON	Sine Wave Triangle Wave ✓ Square Wave Sawtooth Wave (A) ♥0 (Hz) ♥0.5 (A) ♥0	

- 4.11.3. Select Offset, Frequency, and Amplitude.
- 4.11.4. Select if Continue Ramp or Single Ramp.

4.11.5. Set Ramp to ON. The ramp of output will start.

CMUL.		
GMWAssociates Agenet System Closed Loop Control obtware Version: 2.0a. Oct 27, 2007 Exc) STOP INV Setup Log to File Remote Mode Log to File Remote Mode	Operation Menu P Chart Control Chart Control Cha	PID Chart Ramp Mode Table Mode Setup Status
StartLog Merual Log on Merual Log Sample Interval 1 (tec)	Custon:]	Probe Field 0.08317 (T) Sample Field -0.00058 (T)
System Status No error	To Be Add	Output Current 8.31 (A) 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5
Control loop 551 Control loop 525 Error check loop 163		iso -12s -10o -7s -0o -2s -0 -2s -5o -7s -10o 12s 12o

4.12. Table Mode.

At Table mode, software will set the output current step through different set points at time interval defined by user. User can manually input a table or load table from a pre-saved file. The total number of Runs and Time between steps are not part of the table and need to be inputted by user each time the Table Mode was started.

Table Mode can run either on current or field control mode. The letter before the setpoint determines mode, "I" for current control mode and "F" for field control mode.


The table file is a plain txt file. It can be viewed, edited, or generated using an ASCII text editor, such as Notepad. The mode letter (I or F) and setpoint are separated with Tab. The table needs ended by "END", case insensitive. The example table used below has the format as below:

Ι	1
Ι	2
Ι	3
Ι	4
Ι	3
Ι	2
Ι	1
Er	nd

4.12.1. Select from Table Mode tab.

GMW Magnet System Closed Loop Control V2_02.vi Fr File Edit View Project Operate Tools Window Help Image: State Control V2_02.vi Fr Image: State Control V2_02.vi Fr <th>ont Panel - 현고</th> <th></th>	ont Panel - 현고	
GMWAssociates Magnet System Closed Loop Control Software Version: 2.02 Nov 14, 2007 (Esc) STOP HW Setup Log to File Remote Mode Log to File: 8 Remark: Your remark here.	Operation Menu PID Chart Wave Form Table Mode Ramp I	1ode
Start Log Manual Log on Manual Log Sample Interval \$1 (sec) System Status No error	Custom Sample Field, Bp 0.00000 To BE Add Output Current 0.00 To BE Add Output Voltage 0.00	(T)
Loops Read loop 0 Control loop 0 Error check loop 0	-20 -17.5 -15 -12.5 -10 -7.5 -5 -2.5 0 2.5 5 7.5 10 12.5 1	5 17.5 20

4.12.2. Set Table Mode to ON. A Table Mode Setup window will appear.



- Load Table From File. Load the table from a pre-saved file.

- Save Table To File. Save the current table to a user defined file name/directory.


- Start Table Mode. Close the Table Setup window and start table mode.

- Abort Table Mode. Close the Table Setup window and does not start Table mode.

- **Time between steps.** Time span for each step. For field control, this time needs be long enough to allow the field reach the setpoint.

- Total Runs. Total number of runs. Default value is 1, meaning only run through the table value once.

4.12.3. Click on Start Table Mode. The output current will change according to the table.



4.12.4. The table mode will automatically stop after reach the end of table, in the last run. The power supply output will return to the value defined in Operation Menu.


Section 5 MAINTENANCE

EXCITATION CURVES

- Sorensen power supply stability test.
 Kepco power supply stability test.


GMW Associates Electromagnet Stability Plot

Model:	3473-70	Model:	5403EG-50	Engr:	Y.Qin		
SN:	127	Serial Number:	2	Date:	4/13/2007		
Pole Gap:	25mm	Pole Gap:	5mm				
Pole Face:	150mm	Pole Face:	76mm				
Field Meter:	DTM-151(MPT-141-7S), NMR 2025						
Power supply:	Two Kepco BOP 20-20						
Note 1:	DAQ controlling	DAQ controlling two power suppl mag 1 set at 1.2T, mag 2 set at 0.3T(~18A)					
Note 2:	Field mode (closed loop control), ~8 hour test, 1.1T, DTM digital filter ON						
Note 3:	Kepco PS Temp Coefficient: 200ppm/ deg C, 8-hour drift: 200ppm						
Note 4:	DTM Temp Coef	fcient: +/-10ppm/	deg C; Zero +/-3ppm/de	eg C; Cabl	e -3ppm*7m/deg C		



GMW Associates Electromagnet Stability Plot

Model:	3473-70	Model:	5403EG-50	Engr:	Y.Qin		
SN:	127	Serial Number:	2	Date:	4/13/2007		
Pole Gap:	25mm	Pole Gap:	5mm				
Pole Face:	150mm	Pole Face:	76mm				
Field Meter:	DTM-151(MPT-1	DTM-151(MPT-141-7S), NMR 2025					
Power supply:	Two sorensen, 160V/31A						
Note 1:	DAQ controlling	DAQ controlling two power suppl mag 1 set at 1.2T, mag 2 set at 0.25T(~15A)					
Note 2:	Field mode (closed loop control), ~8 hour test, 1.2T, DTM digital filter ON						
Note 3:	Sorensen PS Temp Coefficient: 300ppm/ deg C (0.03%/deg C)						
Note 4:	DTM Temp Coef	fcient: +/-10ppm/	deg C; Zero +/-3ppm/de	g C; Cabl	e -3ppm*7m/deg C		



Section 8 APPENDIX

APPENDIX 1. NATIONAL INSTRUMENTS DAQPad 6015/6016 DATA SHEET APPENDIX 2. NATIONAL INSTRUMENTS USB-DAQ 6251 DATA SHEET APPENDIX 3. SORENSEN SG SERIES POWER SUPPLY DATA SHEET APPENDIX 4. KEPCO BOP SERIES POWER SUPPLY DATA SHEET

Portable High-Performance Multifunction DAQ for USB

NI DAQPad-6015, NI DAQPad-6016

- 16 analog inputs at up to 200 kS/s, 16-bit resolution
- Built-in screw terminal connectors for easier and more
- cost-effective connectivity2 analog outputs for accurate output signals
- 8 or 32 digital I/O lines (5 V TTL/CMOS)

Operating Systems • Windows 2000/XP

Recommended Software

- LabVIEW
- LabWindows/CVI
- Measurement StudioVI Logger

Measurement Services Software (included) • NI-DAQmx

Product	Bus	Analog Inputs ¹	Input Resolution	Sampling Rate	Input Range	Analog Outputs	Output Resolution	Output Rate ²	Output Range	Digital I/O	Counter/Timers	Triggers
NI 6015	USB	16 SE/8 DI	16 bits	200 kS/s	±0.05 to ±10 V	2	16 bits	300 S/s	±10 V	8	2	Digital
NI 6016	USB	16 SE/8 DI	16 bits	200 kS/s	±0.05 to ±10 V	2	16 bits	300 S/s	±10 V	32	2	Digital
¹ SE – singl	1SE – single ended, DI – differential ² System dependent											

Hardware Description

The National Instruments DAQPad-6015 and DAQPad-6016 multifunction data acquisition devices provide plug-and-play connectivity via USB for acquiring, generating, and logging data in a variety of portable and desktop applications. They include built-in screw terminal connectivity so you can easily connect sensors and signals without extra cost. Both devices feature 16-bit accuracy at up to 200 kS/s. The DAQPad-6016 also provides 32 digital I/O lines for applications requiring an extended interface to digital sensors and actuators.

Software Description

NI DAQPad-6015 and DAQPad-6016 devices include NI-DAQmx Measurement Services software with which you can quickly configure and begin taking measurements with your DAQ device. NI-DAQmx provides a seamless interface to LabVIEW, LabWindows/CVI, and Measurement Studio development environments with features such as DAQ Assistant.

Recommended Accessories

The DAQPad-6015 and DAQPad-6016 both have built-in connectivity, so no additional accessories are required.

Common Applications

The DAQPad-6015 and DAQPad-6016 devices are ideal for a number of applications where portability and accurate measurements are essential, such as:

- Portable data logging log environmental or voltage data quickly and easily
- Field monitoring applications
- Embedded OEM applications
- In-vehicle data acquisition
- Academic lab use academic discounts available for quantities of five or more. Visit *ni.com/academic* for details.

Information for OEM Customers

For information on special configurations and pricing, please visit *ni.com/oem*.

Ordering Information

NI DAQPad-6015	779047-01
NI DAQPad-6016	779025-01
Includes NI-DAQmx software	



Portable High-Performance Multifunction DAQ for USB

Specifications-

These specifications are typical at 25 °C unless otherwise stated.

Analog Input

Accuracy Specifications

				Absolute Ac	Relative Accuracy				
	Percent of Reading			Noise + Q	uantization (mV)	Absolute Accuracy	Temperature	Resolut	ion (mV)
Nominal Range (V)	24 hours	1 year	Offset (mV)	Single Point	Averaged	at Full Scale (mV)	Drift (%/°C)	Single Point	Averaged
± 10	0.0658	0.0700	1.8798	0.9330	0.0824	8.984	0.0010	1.0849	0.1085
± 5	0.0158	0.0200	0.9598	0.4665	0.0412	2.003	0.0005	0.5425	0.0542
± 0.5	0.0658	0.0700	0.1158	0.0562	0.0050	0.471	0.0010	0.0663	0.0066
± 0.05	0.0658	0.0700	0.0314	0.0314	0.0031	0.069	0.0010	0.04038	0.0040

Input Characteristics

Number of channels...

Type of ADC ... Resolution.... Maximum sampling rate Input signal ranges (bipolar only)

Device Gain (Software-Selectable)

Device Gain (Software-Selectable)	Kange
0.5	±10 V
1	±5 V
10	±500 mV
100	±50 mV

Input coupling...... DC Maximum working voltage (signal + common mode) ...

Each input should remain within ± 11 V of ground

16 single-ended or 8 differential (software-selectable

per channel)

200 kS/s

Successive approximation 16 bits, 1 in 65,536

Overvoltage protection

Signal Name	Powered Off (V)	Powered On (V)
Al<015>	±15	±25
AI SENSE	±15	±25
FIFO buffer size	4,096 samj	bles
Data transfers		programmed I/O
Configuration memory size	512 words	
Transfer Characteristic	s	
Relative accuracy	±1.5 LSB ty	/pical, ±3.0 LSB maximum
DNL	±0.5 LSB ty	/pical, ±1.0 LSB maximum
No missing codes	16 bits, gu	aranteed
Amplifier Characteristi	cs	
Input impedance		
Normal powered	100 GΩ ir	parallel with 100 pF
Powered off		
Overload		
Input bias current	±200 pA	
Input offset current	±100 pA	
CMRR (DC to 60 Hz)		
Range ± 10 V, ± 5 V		
Range ± 500 mV, ± 50 mV		
Dynamic Characteristic	s	
Small signal (-3 dB) bandwidth	1 425 kHz	
Large signal (17. THD) bandwi	dth 450 kHz	
System noise (LSB _{rms} , includin	g quantization)	
Range	LSB _{rms}	
± 10 V, ± 5 V	0.9	
± 500 mV	1.1	
± 50 mV	6.7	

Settling time to full scale step

Range	Accuracy
± 10 V	±4 LSB, 5µs typical
± 5 V	±2 LSB, 5μs maximum
± 50 mV, ± 500 mV	±2 LSB, 5µs typical

Crosstalk (DC to 100 kHz)

Adjacent channels	-75 dB
All other channels	-90 dB

Ana	log	Output
-----	-----	--------

Accuracy Specifications

	Absolute Accuracy					
	Percent o	f Reading		Absolute Accuracy	Temperature	
Nominal Range (V)	24 hours	1 year	Offset (mV)	at Full Scale (mV)	Drift (%/°C)	
± 10	0.0154	0.0196	1.5680	3.530	0.0005	

ut Characteristics 0

output characteristics	
Number of channels	2 voltage outputs
Resolution	16 bits, 1 in 65,536
Maximum update rate	300 S/s, system dependent
Type of DAC	Double-buffered, multiplying
FIFO buffer size	None
Data transfers	Interrupts, programmed I/O
Transfer Characteristics	
Relative accuracy (INL)	±3 LSB, typical
DNL	±2 LSB, typical

Monotonicity	15 bits
Voltage Output	
Range	±10 V
Output coupling	DC
Output impedance	0.1 Ω , maximum
Current drive	±5 mA, maximum
Protection	Short-circuit to ground
Power-on state	±250 mV

Dynamic Characteristics

Settling time to full-scale step	8 µs to ±1 LSB accuracy
Slew rate	4 V/µs
Noise	360 μV_{rms} , DC to 400 kHz
Glitch energy at mid-scale transition	
Magnitude	±100 mV
Duration	4.0 µs

Digital I/O

Number of channels .. Compatibility

P0.<0..7> Digital logic levels

Level	Minimum (V)	Maximum (V)
Input low voltage	0	0.8
Input high voltage	2.0	5.0
Output low voltage (I _{out} = 2.5 µA)	-	0.4
Output high voltage ($I_{out} = -2.5 \ \mu A$)	4.35	-

8 input/output

5 V TTL/CMOS

Power-on state	Input (high impedance), 1.5 k Ω pull down to D GND
Data transfers	Programmed I/O
Max transfer rate	250 S/s, system dependent

P1.<0..7>, P2.<0..7>, P3.<0..7> (NI 6016 only)

Digital logic levels

Level	Minimum (V)	Maximum (V)
Input low voltage	0	0.8
Input high voltage	2.2	5.0
Output low voltage (I _{out} = 24 mA)	-	0.4
Output high voltage (I _{out} = -13 mA)	3.7	-

Portable High-Performance Multifunction DAQ for USB

Digital Triggers

Purpose	
Analog input	Start, reference, and pause trigger, sample clock
Analog output	Start and pause trigger, sample clock
Counter/timers	Source, gate
Source	PFI <09>
Response	Rising or falling edge, software-selectable
Compatibility Pulse width	5 V TTL 10 ns minimum
Fulse widui	
Calibration	
Recommended warm-up time	15 minutes
Calibration interval	1 year
Onboard calibration reference	
Level	
	actual value stored in EEPROM
Temperature coefficient	
Long-term stability	± 15.0 ppm/sqrt (1000 hours)
Physical	
Enclosure dimensions	17.60 x 14.85 x 3.08 cm
	(8.00 x 6.75 x 1.40 in.)
Weight	
I/O Connectors	Screw terminals
Bus Interface	
Type	LISB 2.0. full speed
19pe	000 2.0, 101 30000
Power Requirement	
9 to 25 V	
Power available at I/O connector	+4.65 to + 5.25 VDC at 1 A
Environment	
Operating temperature	0 to 55 °C
Storage temperature	
Relative humidity	10 to 90%, noncondensing

Certifications and Compliances CE Mark Compliance CE

View additional specifications at ni.com/manuals.

..... 2-wire Power-on state Input (high impedance), 100 k $\!\Omega$ pull up to +5 VDC Data transfers...... Interrupts, programmed I/O Max transfer rate 250 S/s, system dependent Timing I/O Number of channels Up/down counter/timers 2 Frequency scaler..... 1 Resolution

16201011011	
Up/down counter/timers	24 bits
Frequency scaler	4 bits
Compatibility	5 V TTL/CMOS
Digital logic levels	

Level	Minimum (V)	Maximum (V)
Input low voltage	0	0.8
Input high voltage	2.0	5.0
Output low voltage (I _{out} = 5 mA)	-	0.4
Output high voltage (I _{out} = -3.5 mA)	4.35	-

Base clocks available

Handshaking

Up/down counter/timers	20 MHz, 100 kHz
Frequency scaler	10 MHz, 100 kHz
Base clock accuracy	±0.01%
Maximum external source frequency	
Up/down counter/timers	20 MHz
External source selections	PFI <09>
External gate selections	PFI <09>
Minimum source pulse duration	10 ns, edge-detect mode
Minimum gate pulse duration	10 ns, edge-detect mode
Data transfers	
Up/down counter/timer	Interrupts, programmed I/O
Frequency scaler	Programmed I/O

NI Services and Support

NI has the services and support to meet your needs around the globe and through the application life cycle – from planning and development through deployment and ongoing maintenance. We offer services and service levels to meet customer requirements in research, design, validation, and manufacturing. Visit **ni.com/services**.



Training and Certification

NI training is the fastest, most certain route to productivity with our products. NI training can shorten your learning curve, save development time, and reduce maintenance costs over the application life cycle. We schedule instructor-led courses in cities worldwide, or we can hold a course at your facility. We also offer a professional certification program that identifies individuals who have high levels of skill and knowledge on using NI products. Visit ni.com/training.

Professional Services

Our Professional Services Team is comprised of NI applications engineers, NI Consulting Services, and a worldwide NI Alliance Partner Program of more than 600 independent consultants and



integrators. Services range from start-up assistance to turnkey system integration. Visit ni.com/alliance.

OEM Support

We offer design-in consulting and product integration assistance if you want to use our products for OEM applications. For information about special pricing and services for OEM customers, visit ni.com/oem.

Local Sales and Technical Support

In offices worldwide, our staff is local to the country, giving you access to engineers who speak your language. NI delivers industryleading technical support through online knowledge bases, our applications engineers, and access to 14,000 measurement and automation professionals within NI Developer Exchange forums. Find immediate answers to your questions at ni.com/support.

We also offer service programs that provide automatic upgrades to your application development environment and higher levels of technical support. Visit **ni.com/ssp**.

Hardware Services NI Factory Installation Services

NI Factory Installation Services (FIS) is the fastest and easiest way to use your PXI or PXI/SCXI[™] combination systems right out of the box. Trained NI technicians install the software and hardware and configure the system to your specifications. NI extends the standard warranty by one year on hardware components (controllers, chassis, modules) purchased with FIS. To use FIS, simply configure your system online with ni.com/pxiadvisor.

Calibration Services

NI recognizes the need to maintain properly calibrated devices for high-accuracy measurements. We provide manual calibration procedures, services to recalibrate your products, and automated calibration software specifically designed for use by metrology laboratories. Visit **ni.com/calibration**.

Repair and Extended Warranty

NI provides complete repair services for our products. Express repair and advance replacement services are also available. We offer extended warranties to help you meet project life-cycle requirements. Visit ni.com/services.



ni.com • (800) 433-3488

National Instruments • Tel: (512) 683-0100 • Fax: (512) 683-9300 • info@ni.com

High-Performance M Series Multifunction DAQ for USB – 16-Bit, up to 1.25 MS/s, up to 80 Analog Inputs

NI USB-6221, NI USB-6225, NI USB-6229, NI USB-6251, NI USB-6255, NI USB-6259

- Up to 80 analog inputs at 16 bits, 1.25 MS/s (1 MS/s or 750 KS/s scanning)
- Up to 4 analog outputs at 16 bits, 2.8 MS/s (2 µs full-scale settling)
- Up to 48 TTL/CMOS digital I/O lines (up to 32 hardware-timed at up to 1 MHz)
- Two 32-bit, 80 MHz counter/timers
- Analog and digital triggering
- NI-PGIA 2 and NI-MCal calibration technology for improved measurement accuracy
- NI signal streaming for 4 high-speed data streams on USB
- Power supply included
- 1-year warranty
- Additional warranty and calibration services available

Operating Systems

• Windows Vista (32- and 64-bit)/XP/2000

Recommended Software

- LabVIEW
- LabVIEW SignalExpress
- LabWindows™/CVI
- Measurement Studio

Other Compatible Software

- Visual Studio .NET
- C/C++/C#

Measurement Services Software (included)

- NI-DAQmx driver software
- Measurement & Automation Explorer configuration utility
- LabVIEW SignalExpress LE



Family	Bus	Analog Inputs	Resolution (bits)	Max Rate (S/s)	Analog Outputs	Analog Output Resolution (bits)	Max Rate (S/s)	Range (V)	Digital I/O	Clocked DIO ¹	Counter	Resolution (bits)
NI 6251	USB	16	16	1.25 M	2	16	2.86 M	±10	24	8, up to 1 MHz ¹	2	32
NI 6255	USB	80	16	1.25 M	2	16	2.86 M	±10	24	8, up to 1 MHz ¹	2	32
NI 6259	USB	32	16	1.25 M	4	16	2.86 M	±10	48	32, up to 1 MHz ¹	2	32
NI 6221	USB	16	16	250 k	2	16	833 k	±10	24	8, up to 1 MHz ¹	2	32
NI 6225	USB	80	16	250 k	2	16	833 k	±10	24	8, up to 1 MHz ¹	2	32
NI 6229	USB	32	16	250 k	4	16	833 k	±10	48	32, up to 1 MHz ¹	2	32
¹ Correlated	DIO can be o	clocked at up to	1 MHz across the US	SB bus and up to 1	0 MHz using o	nboard regeneration.						

Table 1. Selection Guide for High-Performance M Series Multifunction DAQ for USB

Overview and Applications

With recent bandwidth improvements and new innovations from National Instruments, USB has evolved into a core bus of choice for measurement and automation applications. National Instruments M Series devices for USB deliver high-performance data acquisition in an easy-to-use and portable form factor through USB ports on laptop computers and other portable computing platforms. NI designed a new and innovative patent-pending NI signal streaming technology that enables sustained bidirectional high-speed data streams on USB. The new technology, combined with advanced external synchronization and isolation, helps engineers and scientists achieve high-performance applications on USB.

NI M Series high-performance multifunction data acquisition (DAQ) modules for USB are optimized for superior accuracy at fast sampling rates. They provide an onboard NI-PGIA 2 amplifier designed for fast settling times at high scanning rates, ensuring 16-bit accuracy even when measuring all available channels at maximum speed. All high-performance devices have a minimum of 16 analog inputs, 24 digital I/O lines, seven programmable input ranges, analog and

digital triggering, and two counter/timers. High-speed NI USB-625x and NI USB-622x M Series devices have two-year and one-year calibration intervals, respectively. USB M Series devices are ideal for test, control, and design applications including:

- Portable data logging log environmental or voltage data quickly and easily
- Field-monitoring applications
- Embedded OEM applications
- In-vehicle data acquisition
- Academic lab use academic discounts available

NI Signal Streaming

To optimize the use of the Universal Serial Bus (USB) and deliver high-performance data acquisition, National Instruments created several key technologies to push the limits of USB throughput and latency. NI signal streaming combines three innovative hardware- and softwarelevel design elements to enable sustained high-speed and bidirectional data streams over USB. For more information, visit **ni.com/usb**.



USB M Series for Test

For test, you can use the M Series high-speed analog inputs and 10 MHz digital lines with NI signal conditioning for applications including test, component characterization, and sensor measurement. High-speed USB-625x M Series devices are compatible with the NI SCC signal conditioning platform, providing amplification filtering and power for virtually every type of sensor. This platform is also compliant with IEEE 1451.4 smart transducer electronic data sheet (TEDS) sensors, which offer digital storage for sensor data sheet information.

USB M Series multifunction DAQ devices also complement existing test systems that need additional measurement channels. For higher-channel-count signal conditioning on USB, consider the NI CompactDAQ or SCXI platforms.

USB M Series for Control

USB M Series digital lines can drive 24 mA for relay and actuator control. By clocking the digital lines as fast as 10 MHz (with onboard regeneration), you can use these lines for pulse-width modulation (PWM) to control valves, motors, fans, lamps, and pumps. With four waveform analog outputs, two 80 MHz counter/timers, and four high-speed data streams on USB, M Series devices can execute multiple control loops simultaneously. High-speed USB-625x M Series devices also offer direct support for encoder measurements, protected digital lines, and digital debounce filters. With up to 80 analog inputs, 32 clocked digital lines, and four analog outputs, you can execute multiple control loops with a single device.

You can also create a complete custom motion controller by combining USB M Series devices with the NI SoftMotion Development Module.

USB M Series for Design

For design applications, you can use a wide range of I/O – from 80 analog inputs to 48 digital lines – to measure and verify prototype designs. USB M Series devices and National Instruments LabVIEW SignalExpress interactive measurement software bring benchtop measurements to the PC. With NI LabVIEW SignalExpress, you can quickly create design verification tests. The fast acquisition and generation rates of high-performance USB M Series high-speed devices along with LabVIEW SignalExpress provide fast design analysis. You can convert your tested and verified LabVIEW SignalExpress projects to LabVIEW applications for immediate M Series DAQ use, and bridge the gap between test, control, and design applications.

USB M Series for OEMs

Shorten your time to market by integrating National Instruments OEM products in your design. Board-only versions of USB M Series DAQ devices are available for OEM applications, with competitive quantity pricing and software customization. The NI OEM Elite Program offers free 30-day trial kits for qualified customers. Visit **ni.com/oem** for more information.

Recommended Training and Services

All M Series devices are available with additional warranty and calibration services. For new data acquisition programmers, NI recommends the "Data Acquisition: 7 Steps to Success" tutorial kit. This tutorial kit helps shorten development time for data acquisition applications by describing the various stages of getting started with DAQ including system definition, setup, test, and application programming.

Recommended Software

National Instruments measurement services software, built around NI-DAQmx driver software, includes intuitive application programming interfaces, configuration tools, I/O assistants, and other tools designed to reduce system setup, configuration, and development time. National Instruments recommends using the latest version of NI-DAQmx driver software for application development in National Instruments LabVIEW, LabVIEW SignalExpress, LabWindows/CVI, and Measurement Studio. To obtain the latest version of NI-DAQmx, visit **ni.com/support/daq/versions**. NI measurement services software speeds up your development with features including:

- A guide to create fast and accurate measurements with no programming using DAQ Assistant
- Automatic code generation to create your application in LabVIEW; LabWindows/CVI; LabVIEW SignalExpress; and Visual Studio .NET, C/C++/C#, or Visual Basic using Measurement Studio
- Multithreaded streaming technology for 1,000 times performance improvements
- Automatic timing, triggering, and synchronization routing to make advanced applications easy
- More than 3,000 free software downloads to jump-start your project available at ni.com/zone
- Software configuration of all digital I/O features without hardware switches/jumpers
- Single programming interface for analog input, analog output, digital I/O, and counters on hundreds of multifunction DAQ hardware devices
 M Series devices are compatible with the following versions (or later) of NI application software – LabVIEW, LabWindows/CVI, or Measurement

Studio versions 7.x or LabVIEW SignalExpress 2.x.

Recommended Accessories (Mass-Termination Versions)

Signal conditioning is required for sensor measurements or voltage inputs greater than 10 V. NI SCC products, which are designed to increase the performance and reliability of your data acquisition system, are up to 10 times more accurate than using terminal blocks alone. Refer to Table 2 for more information or visit ni.com/sigcon.

Sensor/Signals (>10 V)		
System Description	Cable	Carrier
SCC Signal Conditioning	SH68-68-EP	SCC
Sensor/Signals (<10 V)		
System Description	Cable	Terminal Block
Screw Terminal (Shielded) ²	SH68-68-EP	SCC-681
BNC Connectivity	SH68-68-EP	BNC-2110
Screw Terminal (Nonshielded) ²	R68-68	SCC-681
Screw Terminal (Shielded) ^{2, 3}	SH68-68-S	SCB-68

¹Includes SCC signal conditioning. ²Consider the integrated screw termination version of the USB DAQ device. ³For use with Connector 1 on USB-6225 and USB-6255 devices.

Table 2. Recommended Accessories

Ordering Information

NI USB-6221	
Screw terminal	779808-0P ¹
NI USB-6225	
Screw terminal	779973-0P ¹
Mass terminal	779974-0P ¹
NI USB-6229	
Screw terminal	779810-0P ¹
NI USB-6251	
Screw terminal	779627-0P ¹
Mass terminal	779694-0P ¹
NI USB-6255	
Screw terminal	779958-0P ¹
Mass terminal	779959-0P ¹
NI USB-6259	
Screw terminal	779628-0P ¹
Mass terminal	779695-0P ¹
Includes NI-DAQmx software.	
¹ P is 1 (U.S. 120 VAC); 2 (Swiss 220 VAC); 3 (Australian 240 VAC)	
4 (Universal Euro 240 VAC); 6 (United Kingdom 240 VAC); 7 (Japa Includes data acquisition driver software, 1 m USB cable, and AC	· · · · · ·
Board-Only Devices for OEM	
	105050.00
NI USB-6221 OEM (qty 1)	
NI USB-6225 OEM (qty 1)	
NI USB-6229 OEM (qty 1)	195959-01
NI USB-6251 OEM (qty 1)	
NI USB-6255 OEM (qty 1)	197201-01

NI USB-6259 OEM (qty 1).....194929-01 Accessories

0-1-1

Cables
SH68-68-EP (shielded)184749-01
SH68-68-S (for USB-62x5)185262-01
R6868 (unshielded ribbon)182482-01
Terminal Blocks and Signal Conditioning Carrier
SCC-2345 carrier777458-01
SCC-68 screw-terminal block for mass termination779475-01
SCB-68 screw-terminal block for mass termination776844-01
BNC-2110 BNC terminal block for mass termination777643-01
Data Acquisition Services
Data Acquisition: 7 Steps to Success779489-01

BUY NOW!

For complete product specifications, pricing, and accessory information, call 800 813 3693 (U.S.) or go to ni.com/usb.

Specifications

>> For complete specifications, see the *NI 622x Specifications* and the NI 625x Specifications manuals at ni.com/manuals.

Specifications listed below are typical at 25 °C unless otherwise noted.

Analog Input

• 1	
Number of channels	
USB-6221/6251	8 differential or 16 single ended
USB-6229/6259	16 differential or 32 single ended
USB-6225/6255	40 differential or 80 single ended
ADC resolution	16 bits
Maximum sampling rate	
USB-6221/6225/6229	250 kS/s single channel,
	250 kS/s multichannel
	(aggregate)
USB-6251/6259	1.25 MS/s single channel,
	1.00 MS/s multichannel
	(aggregate)
USB-6255	1.25 MS/s single channel,
000 0200	750 kS/s multichannel
	(aggregate)
Input coupling	DC
Input coupring	DC
1 5	.10 .5 .1 .0.2.1/
USB-6221/6225/6229	±10, ±5, ±1, ±0.2 V
USB-6251/6255/6259	±10, ±5, ±2, ±1,
	±0.5, ±0.2, ±0.1 V
Maximum working voltage for analog in	
(signal + common mode)	±11 V of AI GND
Input impedance	
Device on	
AI+ to AI GND	>10 G Ω in parallel with 100 pF
AI– to AI GND	>10 G Ω in parallel with 100 pF
Device off	
AI+ to AI GND	820 Ω
AI– to AI GND	820 Ω
Input bias current	±100 pA
Crosstalk (at 100 kHz)	
Adjacent channels	-75 dB
Nonadjacent channels	
USB-6221/6225/6229	-90 dB
USB-6251/6255/6259	-95 dB
Input FIFO size	4,095 samples
Scan list memory	4,095 entries
Data transfers	NI signal streaming on USB,
	programmed I/O
	programmod i O

Analog Triggers (USB-625x Devices Only)

Functions	Start Trigger
	Reference Trigger
	Pause Trigger
	Sample Clock
	Convert Clock
	Sample Clock Timebase

Resolution

Analog Output

Number of channels	
USB-6221/6225/6251/6255	2
USB-6229/6259	4
DAC resolution	16 bits
Maximum update rate	
USB-6221/6225/6229	
1 channel	833 kS/s
2 channels	740 kS/s
3 channels	666 kS/s
4 channels	625 kS/s
USB-6251/6255/6259	
1 channel	2.86 MS/
2 channels	2.00 MS/
3 channels	1.54 MS/
4 channels	1.25 MS/
Timing accuracy	50 ppm o
Timing resolution	50 ns
Output range	
USB-6221/6225/6229	±10 V
USB-6251/6255/6259	±10 V, ±5
	on APFI <
Output coupling	DC
Output impedance	0.2 Ω
Output current drive	±5 mA
Output FIFO size	8,191 san
	among ch
Data transfers	NI signal

Calibration (AI and AO)

Recommended warm-up time	15
Calibration interval	
USB-6221/6225/6229	1 ye
USB-6251/6255/6259	2 ye

Digital I/O/PFI

Static Characteristics

Number of channels USB-6221/6225/6251/6255 USB-6229/6259 Ground reference Direction control Pull-down resistor

Analog edge triggering, analog edge triggering with hysteresis, and analog window triggering 10 bits, 1 in 1,024

s/s per channel S/s per channel S/s per channel IS/s IS/s per channel IS/s per channel 1S/s per channel n of sample rate

±5 V, ±external reference FI <0..1> samples shared channels used NI signal streaming, programmed I/O

minutes

/ear years

24 total, 8 (P0.<07>),
16 (PFI <07>/P1,
PFI <815>/P2)
48 total, 32 (P0.<031>),
16 (PFI <07>/P1,
PFI <815>/P2)
D GND
Each terminal individually
programmable as input or output
50 k Ω typical, 20 k Ω minimum

High-Performance M Series Multifunction DAQ for USB – 16-Bit, up to 1.25 MS/s, up to 80 Analog Inputs

programmed I/O

Waveform Characteristics (Port 0 Only)

Terminals usedUSB-6221/6225/6251/6255USB-6229/6259Port 0 (P0.<0..7>)USB-6229/6259Port/sample sizeUSB-6221/6225/6251/6255Up to 8 bitsUSB-6229/6259Up to 32 bitsWaveform generation (D0) FIFO2,047 samplesDI sample clock frequencyRegenerate from FIFO0 to 1 MHz, system dependentDo to 10 MHzStreaming from memoryData transfersNI signal streaming,

PFI/Port 1/Port 2 Functionality

· · · ·	
Functionality	Static digital input,
	static digital output,
	timing input, timing output
Timing output sources	Many AI, AO, counter, DI,
	DO timing signals
Debounce filter settings	125 ns, 6.425 µs, 2.56 ms,
	disable; high and low transitions;
	selectable per input

General-Purpose Counter/Timers

Number of counter/timers	2
Counter measurements	Edge counting, pulse, semiperiod,
Position measurements	period, two-edge separation X1, X2, X4 quadrature encoding with Channel Z reloading;
	two-pulse encoding
Output applications	Pulse, pulse train with dynamic
	updates, frequency division,
	equivalent time sampling
Internal base clocks	80 MHz
	20 MHz
	0.1 MHz
Base clock accuracy	50 ppm
Inputs	Gate, Source, HW_Arm, Aux, A,
	B, Z, Up_Down
Data transfers	NI signal streaming,
	programmed I/O

Frequency Generator

Number of channels	1
Base clocks	10 MHz, 100 kHz
Divisors	1 to 16
Base clock accuracy	50 ppm
Output can be available on any PFI or R	TSI terminal

Phase-Locked Loop (PLL)

Number of PLLs	1
Reference signal	PXI_STAR,
	PXI_CLK10,
	RTSI <07>
Output of PLL	80 MHz timebase; other signals
	derived from 80 MHz timebase

including 20 MHz and 100 kHz

Any PFI, RTSI, PXI_TRIG,

Software-selectable for

Sample Clock Timebase

Sample Clock Timebase

Gate, Source, HW_Arm, Aux, A,

analog output, digital input, digital output, counter/timer 0,

timebases

PXI STAR

most signals

Start Trigger

Pause Trigger Sample Clock

Convert Clock

Start Trigger Pause Trigger

Sample Clock

B, Z, Up_Down

Sample Clock

Sample Clock

counter/timer 1

V

Reference Trigger

External Digital Triggers

Source
Polarity
Analog input function

Counter/timer functions..... Digital waveform generation (DO) function

Analog output function

Digital waveform	acquisition
(DI) function	

Bus Interface

USB	Hi-Speed USB or full-speed USB
NI signal streaming	4 high-speed data streams;
	can be used for analog input,

Power Requirements

USB power supply requirements 11 to 30 VDC, 20 W

Power Limits

+5 V terminal	1 A max
P0/PFI/P1/P2 and	
+5 V terminals combined	2 A max
Power supply fuse	2 A, 250

Physical Requirements

Enclosure dimensions (includes connect	tors)
Screw termination	26.67 by 17.09 by 4.45 cm
	(10.5 by 6.73 by 1.75 in.)
Mass termination	18.8 by 17.09 by 4.45 cm
	(7.4 by 6.73 by 1.75 in.)

Environmental

Operating temperature 0 to 45 °C Storage temperature...... -20 to 70 °C Humidity..... 10 to 90% RH, noncondensing Maximum altitude...... 2,000 m Pollution degree (indoor use only)..... 2

Safety

This product is designed to meet the requirements of the following standards of safety for electrical equipment for measurement, control, and laboratory use:

- IEC 61010-1, EN 61010-1
- UL 61010-1, CSA 61010-1

Note: For UL and other safety certifications, refer to the product label or visit ni.com/certification, search by model number or product line, and click the appropriate link in the Certification column.

Electromagnetic Compatibility

This product is designed to meet the requirements of the following standards of EMC for electrical equipment for measurement, control, and laboratory use:

- EN 61326 EMC requirements; Minimum Immunity
- EN 55011 Emissions; Group 1, Class A
- CE, C-Tick, ICES, and FCC Part 15 Emissions; Class A

Note: For EMC compliance, operate this device according to product documentation.

CE Compliance

This product meets the essential requirements of applicable European Directives, as amended for CE marking, as follows:

- 2006/95/EC: Low-Voltage Directive (safety)
- 2004/108/EC; Electromagnetic Compatibility Directive (EMC)

Note: Refer to the Declaration of Conformity (DoC) for this product for any additional regulatory compliance information. To obtain the DoC for this product, visit ni.com/certification, search by model number or product line, and click the appropriate link in the Certification column.

Waste Electrical and Electronic Equipment (WEEE)

EU Customers: At the end of their life cycle, all products must be sent to a WEEE recycling center. For more information about WEEE recycling centers and National Instruments WEEE initiatives, visit ni.com/environment/weee.htm.

Environmental Management

NI is committed to designing and manufacturing products in an environmentally responsible manner. NI recognizes that eliminating certain hazardous substances from our products is beneficial not only to the environment but also to NI customers.

For additional environmental information, refer to the NI and the Environment Web page at ni.com/environment. This page contains the environmental regulations and directives with which NI complies, as well as any other environmental information not included in this document.

Management Methods for Controlling Pollution Caused by **Electronic Information Products Regulation (China RoHS)**

电子信息产品污染控制管理办法(中国 RoHS)

🕑 🐵 中国客户 National Instruments 符合中国电子信息产品中限制使用某些有害物质指 令 (RoHS)。关于 National Instruments 中国 RoHS 合规性信息,请登录 ni.com/environment/rohs_china. (For information about China RoHS compliance, go to ni.com/environment/rohs_china.)

NI Services and Support



NI has the services and support to meet your needs around the globe and through the application life cycle – from planning and development through deployment and ongoing maintenance. We offer services and service levels to meet customer requirements in research, design, validation, and manufacturing. Visit **ni.com/services**.

Training and Certification

NI training is the fastest, most certain route to productivity with our products. NI training can shorten your learning curve, save development time, and reduce maintenance costs over the application life cycle. We schedule instructor-led courses in cities worldwide, or we can hold a course at your facility. We also offer a professional certification program that identifies individuals who have high levels of skill and knowledge on using NI products. Visit **ni.com/training**.

Professional Services

Our NI Professional Services team is composed of NI applications and systems engineers and a worldwide National Instruments Alliance Partner program of more than 600 independent consultants and



integrators. Services range from start-up assistance to turnkey system integration. Visit **ni.com/alliance**.

OEM Support

We offer design-in consulting and product integration assistance if you want to use our products for OEM applications. For information about special pricing and services for OEM customers, visit **ni.com/oem**.

Local Sales and Technical Support

In offices worldwide, our staff is local to the country, giving you access to engineers who speak your language. NI delivers industry-leading technical support through online knowledge bases, our applications engineers, and access to 14,000 measurement and automation professionals within NI Developer Exchange forums. Find immediate answers to your guestions at **ni.com/support**.

We also offer service programs that provide automatic upgrades to your application development environment and higher levels of technical support. Visit **ni.com/ssp**.

Hardware Services

NI Factory Installation Services

NI Factory Installation Services (FIS) is the fastest and easiest way to use your PXI or PXI/SCXI combination systems right out of the box. Trained NI technicians install the software and hardware and configure the system to your specifications. NI extends the standard warranty by one year on hardware components (controllers, chassis, modules) purchased with FIS. To use FIS, simply configure your system online with **ni.com/pxiadvisor**.

Calibration Services

NI recognizes the need to maintain properly calibrated devices for high-accuracy measurements. We provide manual calibration procedures, services to recalibrate your products, and automated calibration software specifically designed for use by metrology laboratories. Visit **ni.com/calibration**.

Repair and Extended Warranty

NI provides complete repair services for our products. Express repair and advance replacement services are also available. We offer extended warranties to help you meet project life-cycle requirements. Visit **ni.com/services**.





National Instruments • info@ni.com



© 2007 National Instruments Corporation. All rights reserved. CVI, LabVIEW, Measurement Studio, National Instruments, National Instruments Alliance Partner, NI, ni.com, NI CompactDAQ, NI SoftMotion, RTSI, SCXI, and SignalExpress are trademarks of National Instruments. The mark LabWindows is used under a license from Microsoft Corporation. Other product and company names listed are trademarks or trade names of their respective companies. A National Instruments Alliance Partner is a business entity independent from National Instruments and has no agency, partnership, or joint-venture relationship with National Instruments.



SG Series

Programmable Precision DC Power Supply 5 kW - 150 kW // 40 V - 600 V

FI



ტ



For Further Information or Sales Support Contact Elgar Electronics Corporation www.Elgar.com 858.450.0085

999

Overview

The SG series represents the next generation of high power programmable DC power supplies. Designed for exceptional load transient response, low noise and the highest power density in the industry. The industry leading power density is enhanced by a stylish front air intake allowing supplies to be stacked without any required clearance between units.

At the heart of the SG series is a 5 kW power module. Depending on the output voltage, one to six modules can be configured in a single chassis to deliver 5 kW to 30 kW of power. Combinations of these chassis can then be easily paralleled to achieve power levels up to 150 kW. Paralleled units operate like one single supply providing total system current. Available in two control versions, the SGA has basic analog controls, while the SGI provides intelligent control features

SGI: Advanced Intelligent Control

(Sorensen General purpose Intelligent) The SGI combines onboard intelligent controls with the outstanding power electronics common to all SG family supplies. These controls enable sophisticated sequencing, constant power mode and save/recall of instrument settings. Looping of sequences makes the SGI idea for repetitive testing. An impressive vacuum fluorescent graphical display in eight languages, context sensitive "soft" keys and front panel keyboard simplify programming of the SGI.


SGA: Outstanding Value - Analog Control


(Sorensen General purpose Analog) The SGA, with its industry leading price performance, is available for customers requiring simple front panel analog controls or external control. With the same high performance power electronics as the SGI, the SGA provides essential features like 10- turn potentiometers for setting voltage and current, 3 1/2 digit LED readout plus front panel over-voltage protection (OVP) preview/adjustment and reset.

SGI / SGA Comparison Chart							
Feature	SGA	SGI					
Modular Design	•	•					
Fast Load Transient	•	•					
Parallelable	•	•					
Analog & Digital Summing	Optional	•					
Direct Front Panel V/I Control	•	•					
3½ Digit LED Readout	•						
Graphics Display		•					
Sequencing		•					
Save/Recall Setups		•					
System Power Readouts		•					
Constant Power Mode		•					
IEEE-488.2/RS-232C	Optional	RS-232C Std IEEE-488.2 Optional					
LXI Compliant Ethernet	Optional	Optional					



Compliant Ethernet





Applications

Process Control

Whether you are controlling ion beams for the manufacture of semiconductors, or driving a current through electrolyte for precise control in a plating process, the SG series is an ideal choice with it's small size, reliable modular design and standard analog programmability. Direct control of V and I along with monitoring of the actual voltage and current, provides a simple interface for your PLC or other type of analog controller.

Product Development

Testing & Burn-in of DC-DC converters, laser diodes, automotive and semiconductor components and aircraft flight hardware are just a few of the items being tested using the SG product family. From simple front panel control to complex test sequences for compliance testing, the SG series will keep pace with your changing application needs.

Research

A research environment presents some of the most demanding requirements on your test instrumentation. Equipment that is sufficient today, may not meet the needs of the next project. With the SG series this is no longer a problem. The modular design allows you to easily upgrade to higher power levels in the future, or parallel units to achieve up to 150 kW. With the sophisticated sequencing capability of the SGI model, you can build an infinite variety of test or diagnostic programs and have them execute directly from the power supply.

Automotive Component Test

The 16-bit resolution and Ethernet enabled hardware triggering allows for detailed sequencing associated with battery fluctuation simulation. The tight load regulation capability of the SG series makes it a superior source for validation and acceptance testing and burn-in of automotive components. The 40V models, in particular, provide a full range of testing to simulate battery conditions. Margin testing of components, such as electronic control units (ECU) and electromechanical components, is easily achieved.

Rackmount ATE Systems

The high power density of the SG series makes it ideal for ATE System integration. The wide variety of voltage and current combinations in 3U and 6U heights allows multiple voltage outputs in a small amount of space. The wide variety of control methods possible, allows easy integration into legacy systems as well as high speed systems.

Battery Charging

The SG series provides high accuracy voltage output to optimize battery charging. Battery charging requires high accuracy voltage and stable current output. With the remote interface options, the charging process can easily be automated for volume production.

Applications

Burn-In Materials Research Product Validation

Rackmount ATE Systems

Compliance Testing Process Control Automotive Electronics Battery Charging

See the SFA product brochure for very high current slew rate and low stored energy applications.



Key Features

- **High Power Density:** Up to 15 kW in a 3U / 30 kW in a 6U chassis
- Wide Voltage Range: 0-40V up to 0-600V, in increments of 5 kW from 5 to 30 kW
- Fast Load Transient Response: Protection from undesired voltage excursions
- Low Ripple and Noise: Suitable for the most sensitive applications
- **Parallelable up to 150 kW:** Expandable as your requirement grows
- Modular Design: Upgradeable for the ultimate in investment protection.
- Sequencing: Program custom waveforms
- Easy-to-read Fluorescent Display: SGI supports English, French, German, Italian, Spanish, Chinese, Japanese, and Korean languages
- **16-bit Resolution:** Optional IEEE-488.2 + RS-232C or Ethernet provides precise control
- Ethernet Control: LXI Class C compliant communication with integrated web server
- Direct Relay Control: Control output and sense isolation relays, along with polarity relays. (Ethernet Option Only)
- Hardware Trigger: Ethernet Option Only
- 5 Year Warranty

Product Specifications

-	
Common	
Remote Sense	Load-line loss compensation for models $<= 100$ V is 10% above full scale voltage total (5% per load-line), and models > 100 V is 4% above full scale voltage total (2% per load-line).
Parallel Operation	Up to 5 units may be paralleled for additional cur- rent within the power supply single-unit specifica- tions, with exception of the DC output current set accuracy. Additional paralleled SG units will add 0.3% inaccuracy per unit. To parallel more than 5 units, contact factory.
Series Operation	Up to 2 units (see Output Float Voltage)
Input	
Nominal Voltage 3 phase, 3 wire + ground	208/220 VAC (operating range 187 - 242 VAC) 380/400 VAC (operating range 342 - 440 VAC) 440/480 VAC (operating range 396 - 528 VAC)* *Optional
Frequency	47 – 63Hz
Power Factor	>0.9 typical at 208/220 VAC input >0.78 typical at 380/400 VAC input >0.69 typical at 440/480 VAC input
Protection	$\frac{1}{2}$ cycle ride-though on all three phases, 3 cycle ride through on single phase; missing phase shutdown
Environmental	
Operating Temperature	0 to 50° C
Storage Temperature	-25° C to 65° C
Humidity Range	Relative humidity up to 95% non-condensing, 0° C $-$ 50° C
Altitude	Operating full power available up to 5,000 ft. (~1,500 m), derate 10% of full power for every 1,000 feet higher; non-operating to 40,000 ft. (~12,000 m)
Cooling	Front and side air inlet, rear exhaust. Units may be stacked without spacing.
Regulatory	Certified to UL/CSA 61010 and IEC/EN 61010-1, CE Compliant, Semi-F47 Compliant
Physical	
Dimensions	Width: 19.00" (48.3 cm), Depth 25.0" (63.5 cm) Height: 5-15 kW units: 3U – 5.25" rack mount (13.34 cm) 20-30 kW units: 6U – 10.5" rack mount (26.67 cm)
Weight	3U < 80 lbs. (36 kg) 6U < 160 lbs. (73 kg)
Shipping Weight	See web site for more product & shipping weights.
Programming & Road-h	

Output: Voltage and Current Ranges								
	3U				6U		Ripple 8	& Noise
Power	5 kW	10 kW	15 kW	20 kW	25 kW	30 kW	rms (20 Hz-	р-р (20 Нz-
Voltage			Cur	rent			300 kHz)	20 MHz)
40	125	250	375	500*	625*	750*	20 mV	75 mV
60	83	167	250	333	417	500	20 mV	75 mV
80	63	125	188	250	313	375	20 mV	100 mV
100	50	100	150	200	250	300	20 mV	100 mV
160	31	63	94	125	156	188	25 mV	150 mV
200	25	50	75	100	125	150	25 mV	175 mV
250	20	40	60	80	100	120	30 mV	200 mV
330	15	30	45	61	76	91	30 mV	200 mV
400	12	25	38	50	63	75	30 mV	300 mV
600	8	17	25	33	42	50	60 mV	350 mV

* By way of paralleling 5 kW, 10 kW & 15 kW supplies

Output	
Ripple & Noise (Voltage Mode)	Ripple and noise, typical, measured at full load, nominal AC input. Noise measured with 6 ft. cable, $1\mu f$ at load.
Ripple (Current Mode)	<+/- 0.04% of full scale rms current
DC Voltage Slew Rate	100 ms 5-95% of full scale typical (Contact factory for model specific slew rates)
DC Current Slew Rate	45A / ms typical
Line Regulation	(±10% of nominal AC input, constant load) Voltage Mode: +/- 0.01% of full scale Current Mode: +/- 0.05% of full scale
Load Regulation	(no load to full load, nominal AC input) Voltage Mode: +/- 0.02% of full scale Current Mode: +/- 0.1% of full scale
Load Transient Response	Recovers within 1ms to $+/-0.75\%$ of steadystate output for a 50% to 100% or 100% to 50% load change
Efficiency	87% typical at nominal line and max load
Stability	$\pm 0.05\%$ of set point after 8 hrs. warm-up at fixed line, load and temperature
Temperature Coefficient	0.02%/ C of maximum output voltage rating for volt- age set point 0.03%/ C of maximum output current rating for cur- rent set point
Output Float Voltage	Negative terminal within +/- 150 V of chassis potential.

Programming & Read-back Specifications

	Programming	mming Read-Back / Monitoring			
	Accuracy	Resolution	ion Accuracy Reso		
Front panel Display	SGA: +/- (0.5%fs + 1 digit) SGI, Voltage: +/- 0.1% of full scale SGI, Current: +/- 0.4% of full scale	SGA: 3.5 digits SGI: 4.0 digits	$\begin{array}{l} SGA: +/- \ (0.5\% fs + 1 \ digit)\\ SGI, \ Voltage: +/- \ 0.1\% \ of \ full \ scale\\ SGI, \ Current: +/- \ 0.4\% \ of \ full \ scale \end{array}$	SGA: 3.5 digits SGI: 4.0 digits	Knob control & Display read-back
Remote Analog Interface	Voltage: +/-0.25% of full scale for 0-5 V range, +/-0.5% of full scale for 0-10 V range Current: 0.8% of full scale	NA	+/-1.0% of full scale (0 - 10V)	NA	25-pin D-sub connector (0~5 V or 0~10 V)
Remote Digital Interface	Voltage: +/- 0.1% of full scale, Current: +/- 0.4% of full scale	+/-0.002% of full scale	Voltage: +/- 0.15% of full scale, Current: +/- 0.4% of full scale	+/-0.002% of full scale	RS-232C (Standard on SGI), Optional IEEE- 488.2 and Optional LXI Compliant 10/100 base-T Ethernet (see Options)
OVP	+/- 1% of full scale	+/-0.002% of full scale			Programming range: 5-110% Configured from front panel, remote analog or via optional digital inputs
User I/O	Disconnect & Polarity-reversal relay	Digital 10-pin Molex type connector See www.elgar.com/go/pinouts			
Software	IVI & CVI drivers available under SUPPORT at: www.Elgar.com				

Elgar Electronics Corp. • 9250 Brown Deer Road, San Diego, CA 92121. • 858.450.0085 • 858.458.0267 • Sales@Elgar.com • www.Elgar.com

SERIES **BOP**





BOP, while fully rated power supplies, are also high-powered operational amplifiers with full 4-quadrant, bipolar operation. Their output is capable of both sustained d-c and the replication of arbitrary a-c waveforms.

In Kepco's BOP, the voltage and current outputs can be controlled smoothly and linearly through the entire rated plus and minus ranges, passing smoothly through zero with no polarity switching. BOP are high speed power operational amplifiers that can be used to provide dynamically agile voltage for test and simulation. They are *not* general purpose power supplies. By the nature of the bipolar high speed design, they do not have any energy-storage noise discriminating output capacitors. To realize the full high speed potential of BOP, the load characteristics should be mainly resistive. If the load is capacitive, (>0.1 microfarad), the BOP must be slowed to avoid oscillation. This is accomplished in voltage mode by increasing the value of the feedback capacitance. Special terminals on the user port are available to slow the BOP for optimum current-mode stability into inductive loading.

BOP MODEL TABLE

MODEL(1) (5)	d-c OUTPL E _O max.	JT RANGE I _o max.	CLOSED L VOLTAGE CHANNEL G _V (V/V)	00P GAIN CURRENT CHANNEL G ₁ (A/V)	VOLTAC SERIES R	OUTPUT II Ge mode Series l ⁽²⁾	MPEDANCE Curren Shunt R	NT MODE SHUNT C ⁽³⁾
100 WATT								
BOP 20-5M	± 20V	± 5A	2.0	0.5	80μΩ	20µH	40kΩ	0.05µF
BOP 50-2M	± 50V	± 2A	5.0	0.2	0.5mΩ	100µH	50kΩ	0.05µF
BOP 100-1M	± 100V	± 1A	10.0	0.1	2.0mΩ	200µH	100kΩ	0.05µF
200 WATT								
BOP 20-10M	± 20V	± 10A	2.0	1.0	40μΩ	50µH	20kΩ	0.1µF
BOP 36-6M	± 36V	± 6A	3.6	0.6	120μΩ	50µH	36kΩ	0.1µF
BOP 50-4M	± 50V	± 4A	5.0	0.4	$0.25 m\Omega$	100µH	50kΩ	0.05µF
BOP 72-3M	± 72V	± 3A	7.2	0.3	$0.48 m\Omega$	200µH	72kΩ	0.05µF
BOP 100-2M	± 100V	± 2A	10.0	0.2	1.0mΩ	200µH	100kΩ	0.05µF
BOP 200-1M ⁽⁴⁾	± 200V	± 1A	20.0	0.1	4.0mΩ	1.2mH	200kΩ	0.03µF
400 WATT								
BOP 20-20M	± 20V	± 20A	2.0	2.0	20μΩ	50µH	20kΩ	0.2µF
BOP 36-12M	± 36V	± 12A	3.6	1.2	60μΩ	50µH	36kΩ	0.2µF
BOP 50-8M	± 50V	± 8A	5.0	0.8	125μΩ	100µH	50kΩ	0.15µF
BOP 72-6M	± 72V	± 6A	7.2	0.6	240μΩ	200µH	72kΩ	0.1µF
BOP 100-4M	± 100V	± 4A	10.0	0.4	500μΩ	200µH	100kΩ	0.1µF

(1) For factory installed digital interfaces add appropriate suffix. See page 55.

(2) For determining dynamic impedance in voltage mode.

(3) For determining dynamic impedance in current mode.

(4) Same size as 400W models.

(5) To specify digital display, substitute the suffix letter "D" for the suffix letter "M."

FEATURES

- Source and sink 100% of their current rating. See Figure 1.
- Separate control circuits for voltage and current with automatic crossover to current and voltage limits.
- All controls and flag signals accessible through a 50-terminal userport at the rear.
- Zeroable preamplifier available for scaling and summing external signals.
- Optional digital displays. Specify by substituting the suffix "D" in place of the "M."



For high power bipolar power supplies, see Series BOP High Power, page 44.

For high voltage bipolar power supplies, see Series BOP-HV, page 56.

The tabulation of the effective series resistance and inductance in voltage mode, and the effective shunt resistance and shunt capacitance in current mode, is done to allow a calculation of the output impedance versus frequency.

BOP accept plug-in cards for remote digital control

- BIT 4882 provides 12-bit IEEE 488.2 talklisten control with SCPI support.
- BIT 4886 provides 16-bit IEEE 488.2 talklisten control with SCPI support.
- BIT TMA-27 connect BOP to Kepco's single-address multiple instrument serial bus for long range (>300m) control from IEEE 488.2, RS 232 or VXI-based hosts.
- BIT 488B or BIT 488D offer listen-only GPIB support in binary or Hex format.

Cards may be factory installed. See page 55 for appropriate suffix designations.





KEPCO, INC. • 131-38 Sanford Avenue • Flushing, NY 11352 USA • Tel: (718) 461-7000 • Fax: (718) 767-1102 Email: hq@kepcopower.com • www.kepcopower.com/bop.htm







TABLE 1		
Source current measured worst case,	125V	a-c.

MODEL	CURRENT (Amps)			
BOP 20-5M	2.6			
BOP 20-10M	5.5			
BOP 20-20M	11.0			
BOP 36-6M	5.1			
BOP 36-12M	10.6			
BOP 50-2M	2.6			
BOP 50-4M	4.8			
BOP 50-8M	9.5			
BOP 72-3M	5.0			
BOP 72-6M	10.8			
BOP 100-1M	2.6			
BOP 100-2M	4.8			
BOP 100-4M	9.2			
BOP 200-1M	5.5			





BOP GENERAL SPECIFICATIONS

.

BOP GENERAL SP	ECIFICATIONS				
SPECIFICATION	RATING/DESCRIPTION	CONDITION			
a-c Voltage	95-113, 105-125, 190-226, 210-250V a-c	User selectable			
Current	See Table 1	Max load, 115V a-c			
Frequency	47-65Hz	Range			
OUTPUT					
d-c Output	Bi-direction, series pass	Transistor ⁽¹⁾			
Type of stabilizer	Automatic crossover	Voltage/current			
Voltage	0 to 100% of rating (bipolar)	Adjustment range			
Current	0 to 100% of rating (bipolar)	for temp 0-55°C			
Sink	See source/sink plot	Duty cycle			
Error Sense	0.5V per load wire	Voltage allowance			
Isolation Voltage	500V d-c or peak	Output to ground			
Leakage Current	<5 microamperes	rms at 115V a-c 60Hz			
Output to Ground	<50 microamperes	p-p at 115V a-c 60Hz			
Series Connection	500V	Max voltage off ground			
Parallel Connection	Current sharing	Use master-slave connection			
OVP	Not available				
CONTROL					
Type Voltage	Variable input, fixed gain				
Current Voltage/ Current	1,,				
Local	10-turn zero-center pot				
Remote Analog	-10V to +10V				
Local Digital	Serial bus or GPIB or VXI	Optional internal BIT card			
Remote Digital	Use SN or SNR interface				
Bounding	±Volt/current local	Four screwdriver trimmers			
0	±Volt/current remote	0 to 10 volts			
Dynamics	See dynamic spec table	Fast only			
User Amplifiers	Uncommitted gain 20K	Two provided			
References	±10 volts, 1mA	Two provided			
Options (built-in)	GPIB hex card	Suffix -488B			
For user added	GPIB BCD card	Suffix -488D			
card refer to "BIT"	Long range serial card/VXI	Suffix -TMA			
models page 55	Talk-listen 4882 card (SCPI)	Suffix -4882			
	Talk-listen 4886 card (SCPI)	Suffix -4886			
	Serial RS 232	Suffix -232			
MECHANICAL					
Input Connection	Detachable IEC type 3-wire	All models			
Output	Front signal/output	Binding posts			
Connections	Rear user port	50-terminal connector			
	Rear output	Barrier strip			
Meters	Two 2½ ["] horiz., 2% zero center analog	Front panel			
Indicators	Four LEDs	Voltage/Current/Bounding			
Mounting	Use RA 37 rack adapter	³ ⁄ ₄ rack size			
(in std 19" racks)	Mounting "ears" supplied	Full rack size			
Cooling	Forced air	Exhaust to rear			
Dimensions inches (HxWxD) add mm	5 ⁷ ⁄ ₃₂ x 12 ¹ 7⁄ ₃₂ x 17% ₄ 132.6 x 318.3 x 435.4	³ ⁄ ₄ rack size			
21/2" to rear inches for connector mm	5 ⁷ / ₃₂ x 19 x 20 ⁵ / ₆₄ 132.6 x 482.6 x 510	Full rack size			
Finish; Fed Std 595	Light gray, color 26440	Front panel			
Weight	47lb (21.4Kg)	³ / ₄ rack size (100W)			
(packed for	53lb (24.1Kg)	³ / ₄ rack size (100W)			
shipment)	76lb (34.5Kg)	Full rack size			

(1) 200V model uses FET.

BOP DYNAMIC SPECIFICATIONS										
MODEL	BANDWIDTH (d-c to f _{-3dB}) KHz (minimum) Mode V I		RISE & FALL TIME 10%-90% µsec (maximum) Mode V I I		LARGE SIGNAL FREQUENCY (min) RESPONSE, KHz Mode V I		SLEWING RATE (minimum) Mode V I I		RECOVERY STEP LOAD µsec (maximum) Mode V I	
100 WATT										
BOP 20-5M	18	12	20	30	17	13	5V/µsec	0.15A/µsec	25	10
BOP 50-2M	18	12	20	30	17	13	5V/µsec	0.15A/µsec	25	10
BOP 100-1M	18	11	17	22	18	11	11V/µsec	70mA/µsec	40	25
200 WATT	•							•	•	
BOP 20-10M	18	6	20	60	17	7	2V/µsec	0.4A/µsec	80	20
BOP 36-6M	16	13	20	27	15	14	3V/µsec	0.5A/µsec	50	35
BOP 50-4M	23	14	14	25	15	11	4.5V/µsec	0.25A/µsec	40	30
BOP 72-3M	20	15	18	26	17	12	10V/µsec	0.15A/µsec	30	30
BOP 100-2M	22	15	18	26	17	12	10V/µsec	0.15A/µsec	30	30
BOP 200-1M	4.0	2.5	110	150	4.0	2.5	5V/µsec	15mA/µsec	150	120
400 WATT										
BOP 20-20M	9.5	10	35	35	8	10	1V/µsec	1.25A/µsec	100	75
BOP 36-12M	20	10	16	30	19	10	4V/µsec	0.75A/µsec	50	30
BOP 50-8M	24	10	14	35	24	11	7.5V/µsec	0.5A/µsec	40	30
BOP 72-6M	19	9.5	18	40	20	11	9V/µsec	0.4A/µsec	50	20
BOP 100-4M	18	14	22	30	16	10	10V/µsec	0.25A/µsec	40	30

BOP STATIC SPECIFICATIONS

INFLUENCE QUANTI	TY	VOLTAGE TYPICAL		EFFECTS(1) CURREN TYPICAL	r mode Maximum	PREAMI OFFS ΔE _{io}	REFERENCE ± 10V	
Source (minma	Source (minmax.) <0		0.001%	<0.002%	0.005%	<5µV	<1nA	<0.0005%
Load (NL-FL)	Load (NL-FL) <0.00		0.002%	<0.5mA	1mA	—	—	<0.0005%
Time (8-hour drift)		<0.005%	0.01%	<0.01%	0.02%	<20µV	<1nA	<0.005%
Temp., per °C		<0.005%	0.01%	<0.01%	0.02%	<20µV	<1nA	<0.005%
Ripple and Noise (2)	rms	<1mV	3mV ⁽⁵⁾	<0.01%	0.03%	—	_	—
	p-p ⁽³⁾	<10mV	30mV ⁽⁵⁾	<0.1%	0.3%	_	_	

(1) Output effects, expressed as a percentage, are referred to the maximum rated output voltage or current.

(2) Measured with the common terminal grounded so that the common mode current does not flow through the load.

(3) Peak-to-peak ripple is measured over a 20Hz to 10MHz bandwidth. (4) The output effect can be calculated by the relationship: $\Delta E_{0} = \pm \Delta E_{r} (R_{f}/R_{i}) \pm \Delta E_{i0} (1+R_{f}/R_{i}) \pm \Delta I_{i0} (R_{f}) \text{ where } R_{f} \text{ is the feedback resistor,}$ and R_{i} is the input resistor from the reference, E_{r} . (5) For BOP 200-1M the maximum ripple and noise is 5mV rms and 50mV p-p.


The tabulated offsets, more particularly their change as a function of source, time and temperature, allow a user to calculate performance of the uncommitted amplifier(s) with user specified input and feedback components. The formula for this is given in the static specifications table footnote.



KEPCO, INC. • 131-38 Sanford Avenue • Flushing, NY 11352 USA • Tel: (718) 461-7000 • Fax: (718) 767-1102 Email: hg@kepcopower.com • www.kepcopower.com/bop.htm

OUTLINE DIMENSIONAL DRAWINGS

Fractional dimensions in light face type are in inches, dimensions in bold face type are in millimeters. Tolerance: \pm 1/64" (0.4) between mounting holes \pm 1/32" (0.8) other dimensions



BOP 3/4 rack width power supplies can be rack mounted using RA 37. See page 77.

KEPCO, INC. • 131-38 Sanford Avenue • Flushing, NY 11352 USA • Tel: (718) 461-7000 • Fax: (718) 767-1102 Email: hq@kepcopower.com • www.kepcopower.com/bop.htm