


---

# LM solenoid magnet systems

---



Author: Dr. Taotao Huang

Document Ref: LM2015v2

Email: [t.huang@hts-110.com](mailto:t.huang@hts-110.com)

Initial Issue Date: 15/9/2015

## Content

|       |   |    |
|-------|---|----|
| 1     | System Description .....                    | 3  |
| 1.1   | Magnet.....                                 | 3  |
| 1.2   | Specification .....                         | 3  |
| 1.3   | Magnet Layout & Performance.....            | 4  |
| 1.3.1 | Magnetic Field Profile .....                | 5  |
| 1.3.2 | Fringe Field.....                           | 9  |
| 1.4   | Power Supply .....                          | 14 |
| 1.5   | Cryocooler .....                            | 14 |
| 1.6   | Magnet Field Control Software.....          | 14 |
| 1.7   | Magnet System Supervisor .....              | 15 |
| 1.8   | Yoke .....                                  | 15 |
| 2     | Infrastructure requirements (expected)..... | 16 |
| 3     | Pre-delivery Tests .....                    | 17 |
| 4     | Technical Information .....                 | 18 |
| 4.1   | Magnet System Supervisor (MSS) .....        | 18 |
| 4.2   | Key MSS Interfaces .....                    | 18 |
| 4.3   | Power Supply Data Sheet .....               | 19 |
| 4.4   | Cryocooler Data Sheet .....                 | 20 |
| 4.5   | Key Magnet Interfaces.....                  | 21 |
| 4.6   | Outline drawing (LM-53-2T) .....            | 22 |

### **Please Note**

*All designs and/or information embodied in this document are copyright © of HTS-110 and shall not be provided to third parties. Designs and information may be based on techniques protected by patents owned, be subject of patent applications by HTS-110, or be subject to other forms of Intellectual Property protection. Use of these designs or information by the recipient (or other parties associated or not associated with the recipient) is subject to agreement of licensing and/or other commercial arrangements between said parties and HTS-110. Where fees have been paid (or other considerations made) by recipient for a design study, preparation of samples, testing of samples or other work by HTS-110 personnel, acceptance of such payments by HTS-110 does not confer title or any rights to use designs and/or information embodied herein or otherwise supplied by HTS-110.*

# 1 System Description

## 1.1 Magnet

HTS-110's LM series will provide a maximum field of 1-3T at the centre of the Ø40-80mm room temperature bores. LM series are cryogen-free and are cooled by means of a cryocooler with an associated compressor.


LM series will be provided with a bipolar Kepco power supply which allows the field to be set at any desired level or ramped continuously as required. The power supply is provided with interfaces which allow it to be operated via digital communications from a PC.

LM series are provided with monitoring electronics ("Magnet System Supervisor") to ensure safe, reliable operation.

## 1.2 Specification

| Model                     | Bore Size<br>( mm ) | Peak Field<br>( T ) | Height<br>( mm ) | Cryostat OD<br>( mm ) | Uniformity, 5mm DSV<br>( rms ) ( % ) |
|---------------------------|---------------------|---------------------|------------------|-----------------------|--------------------------------------|
| LM-40-1T                  | 40                  | 1.4                 | 45               | 254                   | 0.4                                  |
| LM-40-2T                  | 40                  | 2                   | 55               | 254                   | 0.3                                  |
| LM-40-3T                  | 40                  | 3                   | 100              | 254                   | 0.2                                  |
| LM-53-1T                  | 53                  | 1.1                 | 45               | 254                   | 0.3                                  |
| LM-53-2T                  | 53                  | 2                   | 65               | 254                   | 0.2                                  |
| LM-53-3T                  | 53                  | 3                   | 120              | 254                   | 0.15                                 |
| LM-80-1T                  | 80                  | 1                   | 55               | 254                   | 0.2                                  |
| LM-80-2T                  | 80                  | 2                   | 100              | 254                   | 0.15                                 |
| LM-80-3T                  | 80                  | 3                   | 160              | 254                   | 0.05                                 |
| Maximum operating current |                     |                     |                  | 125A                  |                                      |

### 1.3 Magnet Layout & Performance



The magnets can be oriented vertically and horizontally. Please note that support brackets are optional and not included in this supply.

### 1.3.1 Magnetic Field Profile




Figure 1 Field profile of LM-40-1T in the axial direction







Figure 2 Field profile of LM-40-2T in the axial direction




**Figure 3 Field profile of LM-40-3T in the axial direction**




**Figure 4 Field profile of LM-53-1T in the axial direction**




**Figure 5 Field profile of LM-53-2T in the axial direction**



**Figure 6 Field profile of LM-53-3T in the axial direction**




**Figure 7 Field profile of LM-80-1T in the axial direction**



**Figure 8 Field profile of LM-80-2T in the axial direction**






**Figure 9 Field profile of LM-80-3T in the axial direction**

Figure 1-9 shows the fields in the axial direction.

### 1.3.2 Fringe Field



**Figure 10 Fringe field of LM-40-1T**




Figure 11 Fringe field of LM-40-2T




Figure 12 Fringe field of LM-40-3T




Figure 13 Fringe field of LM-53-1T




Figure 14 Fringe field of LM-53-2T




Figure 15 Fringe field of LM-53-3T






Figure 16 Fringe field of LM-80-1T



**Figure 17 Fringe field of LM-80-2T**



**Figure 18 Fringe field of LM-80-3T**

50 Gauss (red) and 5 Gauss (blue) lines of the magnets are shown in Figure 10-18.

## **1.4 Power Supply**

KEPCO BOP power supplies are normally used with HTS-110 magnets. These are 4-quadrant supplies offering smooth sweeps through a zero field (no switching). The power supplies have communication interfaces for remote operation and are connected to the Magnet System Supervisor for safety shut-down. The magnets are driven by one KEPCO BOP 6-125MG power supply.

## **1.5 Cryocooler**

The magnet coils are kept at operating temperature using a cryocooler and associated compressor. Specifications of this cooler are included in the Technical Information Section.

## **1.6 Magnet Field Control Software**

The relationship between field and current is determined when the magnets are tested at HTS-110. The power supply can be controlled from a PC on which the client can run its own programs to control the magnet field.

## **1.7 Magnet System Supervisor**

All HTS-110 magnets are supplied with an in-house magnet system supervisor (MSS). The MSS monitors temperature at several points within each coil pack as well as voltage across the coils. Given the high thermal stability of HTS magnets, abnormal operating condition can be readily detected and acted on before a quench occurs. If a pre-quench condition is detected a shut-down command is sent to the power supply and the stored energy dumped through a resistor based dump circuit.

The dump circuit is sized to ensure complete energy dump within the required time from the sequence trigger. Approximate decay time in case of dump is in the range 1-5 seconds.

## **1.8 Yoke**

An aluminium cryostat integrated with returning yoke forms the mechanical chassis for the system. This outer yoke is manufactured from low carbon steel (AISI 1006) with a nickel coating for corrosion resistance.

## 2 Infrastructure requirements (expected)

For installation an oil-free vacuum pump (not supplied by HTS-110) is required. It is also recommended that a vacuum pump is available on site so that the magnet vacuum can be reconditioned from time to time.

|                          |  |
|--------------------------|--|
| Cryocooler & Compressor  | 1 phase 200V, 220V,240V, 50Hz or<br>200V, 208V, 220V 60Hz, 11A (Max) |
| Cooling Water            | 5 litre per minute @ 18°C  |
| Magnet System Supervisor | 100-240 VAC,47/63HZ, 1ph 0.25A (Max)                                 |
| Magnet Power Supply      | 176-264 VAC, 47/63 Hz, 1ph, 9.5 A (Max)                              |
| Current Terminal Heaters | 110-230 V : 50/60 Hz, 1ph 0.1A (Max)                                 |



### **3 Pre-delivery Tests**

1. Vacuum integrity
2. Cool-down time
3. Static monitor tests
  - a. Disable signal correctly switches off the power supply
  - b. Power-supply disabled when voltage limits exceeded
  - c. Power-supply disabled when temperature limits exceeded
  - d. Power supply disabled if magnet monitor cable is removed/monitor switched off
4. High current monitor/dump-circuit test - Power supply disabled and magnet safely ramps down
5. Fringe field - Radial and axial 5 gauss line
6. B vs I excitation curves for magnet
7. Field maps if required
8. Single non-linear cycle at fastest stable rate.
  - Cycle time and maximum temperature excursion recorded.

## 4 Technical Information

### 4.1 Magnet System Supervisor (MSS)

The Magnetic System Supervisor (MSS) is designed to monitor the conditions of the HTS magnet and automatically disable the power supply if the magnet begins to operate abnormally. The MSS can control the power supply attached to the magnet system via the auxiliary communication port using a standard RS-232 cable to create a fully autonomous system.


The MSS fits inside 19" Rack, it is 1U tall, 10" deep and weighs 3kg. The MSS requires 47/63Hz 100-240V, 0.4A max, single-phase power.

#### Features

- Five dedicated temperature channels: for simultaneously measuring up to five cryogenic temperature sensors.
- Two dedicated voltage channels: for simultaneously measuring up to two fully differential voltages.
- PC communication output: for remote control/monitoring of the system via a computer.
- Interlock relay output: for disabling the power supply in the event of a quench.
- Aux relay output: for switched control output. (optional)
- Status lamp: a multi-colour led for visual status update.



### 4.2 Key MSS Interfaces



### 4.3 Power Supply Data Sheet



Detailed input characteristics:


| SPECIFICATIONS          |              | RATING/DESCRIPTION                      | CONDITION  |
|-------------------------|--------------|---|--|
| a-c Voltage             | nominal      | 230V a-c                                | Single phase   |
|                         | range        | 176-264V a-c                            |  |
| Frequency               | nominal      | 50/60 Hz                                |  |
|                         | range        | 47-63 Hz                                |  |
| Current                 | 176V a-c     | 9.5A (7.5A)*                            | Maximum  |
|                         | 264V a-c     | 6.4A (4.4A)*                            | Maximum  |
| Power Factor            | source       | 0.99                                    | Nominal output power   |
|                         | sink         | 0.97                                    |  |
| Efficiency              |              | 65% (56%)*                              | Minimum when sourcing  |
| Switching Frequency     |              | 70 KHz $\pm 5\%$<br>(50KHz $\pm 5\%$ )* | Active PFC for<br>source and sink  |
| EMC Compliance          |              | EN61326-1 (1997)                        | Class A equipment  |
| EMC Immunity            | ESD          | EN61000-4-2                             | Electrostatic discharge  |
|                         | Radiated RF  | EN61000-4-3                             |  |
|                         | EFT          | EN61000-4-4                             | Electrical fast transient/burst  |
|                         | Surges       | EN61000-4-5                             |  |
|                         | Conducted RF | EN61000-4-6                             |  |
| EMC Emissions           | Conducted    | EN61000-3-2                             | Harmonics  |
|                         |              | EN61000-3-3                             | Fluctuation and flicker  |
|                         | Conducted    | EN55011/CISPR11                         | 0.15 to 30 MHz   |
|                         | Radiated     | EN55011/CISPR11                         | 30 to 1000 MHz   |
| Leakage Current         |              | 3.5 mA                                  | 230V a-c, 47-63 Hz   |
| Insulation Coordination | Input        | Installation Category II                | For TN or TT power system  |
|                         |              | Overvoltage Category II                 |  |
|                         | Output       | Installation Category II                | Maximum 300V common mode<br>voltage between output terminals<br>and chassis ground |
|                         |              | Overvoltage Category II                 |  |
| Pollution Degree        |              | 2                                       |  |

## 4.4 Cryocooler Data Sheet

### Coolstar 0/40 cold head

|  | Coolstar 2/9      | Coolstar 6/30     | Coolstar 0/12     | Coolstar 0/40     |
|--|-------------------|-------------------|-------------------|-------------------|
| <b>Base Temperature (Unloaded)</b>                           | 10 K              | 10 K              | 30 K              | 30 K              |
| <b>Nominal Cooling Power at 20 K</b>                         |                   |                   | -                 | -                 |
| Normal Speed (72 rpm)  | 2 W               | 6 W               | -                 | -                 |
| High Speed (90 rpm)  | 2.4 W             | 7 W               | -                 | -                 |
| <b>Nominal Cooling Power at 77 K</b>                         |                   |                   |                   |                   |
| Normal Speed (72 rpm)  | 9 W               | 28 W              | 12 W              | 42 W              |
| High Speed (90 rpm)  | 11 W              | 30 W              | 17 W              | 57 W              |
| <b>Time to cooldown bare coldhead</b>                        |                   |                   |                   |                   |
| to 20 K  | <16 min           | <18 min           |                   |                   |
| to 77 K  |                   |                   | <11 min           | <11 min           |
| <b>Increase in cooldown time for each 100 g Copper added</b> | 14 min            | 5.2 min           | 6 min             | 2 min             |
| <b>Weight</b>  | 2.5 kg            | 6.1 kg            | 2.4 kg            | 5.6 kg            |
| <b>Dimensions</b>  | 278 x 129 x 150mm | 344 x 155 x 207mm | 185 x 123 x 150mm | 219 x 155 x 207mm |



#### 0/40 COLDHEAD



### Cryodrive 2.0 compressor

|   |   |
|---|---|
| <b>Electrical Power (<math>\pm 6\%</math>)</b>        | Only single-phase power required<br>200 V, 220 V, 240 V 50 Hz or 200 V, 208 V, 220 V, 60 Hz |
| <b>Typical Running current:</b>                       |   |
| Cryodrive 2.0   | 11 A  |
| Cryodrive 3.0   | 13 A  |
| <b>Typical Helium Pressure</b>                        |   |
| Running   | 22.5 Bar  |
| Static  | 16.5 Bar  |
| <b>Operating Temperature</b>                          | 4-38°C  |
| <b>Typical Cooling Water Flow</b>                     | 5 l/min at 18°C   |
| <b>Weight + Dimensions</b>                            | 80 kg, 451 x 556 x 448mm  |
| <b>Standard Helium Hose and Coldhead Cable Length</b> | 3 m   |
| <b>Cooling Water Connection</b>                       | 1/2" ID nozzle  |

## 4.5 Key Magnet Interfaces



#### 4.6 Outline drawing (LM-53-2T)

