

OVERVIEW

The 3480 dipole electromagnet is a light weight versatile system that can provide fields approaching 4Tesla. At 34kg this magnet can easily be moved between applications and can be operated in any orientation.

The system is intended to be compatible with systems that already utilize GMW's 3470 electromagnet.

The 3480 is shipped with a standard set of poles that optimize maximum field, but several pole options are available. GMW can also design custom poles that achieve a specific performance. Poles are interchangeable and are available with an axial access bore.

3480 Electromagnet

Features

- Weight: 34kg
- Peak fields to 3.6 T (5mm pole face diameter at 2mm gap)
- Any Mounting Orientation
- Fast Cycle Times

Applications

- Hall Effect Studies
- Magneto-Optical Studies
- Laboratory Experiments
- Magnetic Separation Experiments
- Magnetic Alignment of Small Components
- Chemical Reaction Rate Studies

GMWAssociates

- 🕀 www.gmw.com
- 🔀 sales@gmw.com
- **\$** +1-650-802-8292
- 955 Industrial Road
 San Carlos, California, USA

Model 3480 General Specifications

Mechanical

Dimensions	285mm W x 197mm Dx 213mm H			
Weight (excluding hoses and water)	34kg			
Standard Pole Face Diameter	5mm, 16mm, 20mm, 40mm			
Pole Gap	0 to 75mm			
Coil Spacing	47 to 75mm			

Coils (series connected)

Resistance (20°C)	1.46Ω
Max. Resistance	1.76Ω
Low Current Inductance	0.23H
High Current Inductance	0.13H
Anticipated max. Sinusoidal frequency(1T)	10Hz
Max. Continuous Power (water)	35A, 60V
Max. Continuous Power (air)	7A, 12V
Peak Power (Sinusoid) (water)	50A, 83V
Peak Power (Sinusoid) (air)	10A, 17V
Peak Power (Triangle Wave) (water)	60A, 100V
Peak Power (Triangle Wave) (air)	12A, 20V
Water Cooling (supply 18°C)	4liters/min, 1.0bar (1US GPM, 15psid)

Safety

Overtemperature Interlock	Selco 802L-065 thermostat, mounted onto each cooling plate, wired in series. Contacts below 65°C
Diameter Sphere Containing 5G-surface ("fringe field")	1000mm





Figure 1: Comparison of performance of different pole pairs. The data reflects performance at 35A with coils at maximum spacing.











Figure 2: Current versus voltage curve for all windings connected in series. The departure from linearity is due to coil heating.





Figure 3: Peak achievable field against available power for standard pole (16mm pole face diameter) at 10mm pole gap and minimum coil spacing. This plot assists the selection of power source for desired performance.





Application Note: Ferromagnetic Resonance

A custom pole is designed to create peak field across a 2mm pole gap for the purposes of FMR analyses. A comparison with the standard 16mm pole is shown below.



Figure 4: Field versus pole gap for two pole options. Excitation is 35A.

For any pole gap an optimized pole is available, on request, which will provide best performance for any specific pole gap.



Field Uniformity in the Median Plane



Figure 5: Field profile in the median plane for the 5mm pole face diameter pole. Excitation current is 35A and coils are at minimum spacing. The field uniformity is better than $\pm 0.5\%$ over Ø3mm.

- Peak continuous field 3.71 Tesla
- Peak pulsed field 3.85 Tesla
- Field uniformity <±0.5% over Ø3mm
- Higher uniformity field available on request.



Figure 6: Performance of optimized pole for 2mm pole gap. Continuous operation is possible up to 35A with a triangle waveform possible up to 60A peak.



Application Note: EPR Spectroscopy

GMW offers poles optimized for EPR applications. Two examples are shown below for pole face spacing of 14mm. Each pole give good uniformity across the excitation range for frequencies up to 40GHz but each is optimized to give excellent field uniformity at either K-band or the Q-band. Similar custom poles are available on request for any frequency in the 0-40GHz frequency band.

Waveband	L	S	С	Х	Р	К	Q	U	V
λ [mm]	300	100	75	30	20	12.5	8.5	6	4.6
f [GHz]	0.8-1.2	3.4-3.8	4	9-10	15	24	34	50	65
B [T]	0.035	0.13	0.14	0.34	0.54	0.85	1.22	1.8	2.3



